**cercon**  
smart ceramics

## Cercon<sup>®</sup> eye

Gebrauchsanweisung  
Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso  
Instruções de utilização  
Gebruiksaanwijzing

cercon® eye

DeguDent

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The text and images for this operating manual were prepared with great care. However, typographical or other incorrect specifications may occur. Please observe that DeguDent GmbH accepts no liability for typographical or specification errors.

<b>1</b>	<b>Indications for use</b> .....	<b>20</b>
1.1	Indications for use .....	20
1.2	About this instruction for use .....	20
1.3	Organisational elements used in this manual .....	20
1.4	Device type and year of manufacture .....	20
1.5	Manufacturer, Distributor .....	20
1.6	Copyrights and property rights .....	20
<b>2</b>	<b>Contraindications</b> .....	<b>21</b>
2.1	Contraindications .....	21
2.2	Transport safety information .....	21
2.3	Operating safety information .....	22
<b>3</b>	<b>Warnings</b> .....	<b>22</b>
3.1	Warnings .....	22
<b>4</b>	<b>Precautions</b> .....	<b>23</b>
4.1	Precautions .....	23
4.2	Front side .....	23
4.3	Rear side .....	23
4.4	Interior with rotary plate .....	24
4.5	Modelling table with gauge marks .....	24
<b>5</b>	<b>Advers reactions</b> .....	<b>25</b>
5.1	Advers reactions .....	25
5.2	Setting up .....	25
5.3	Electrical connection .....	25
<b>6</b>	<b>Scope of delivery</b> .....	<b>26</b>
6.1	Including in the delivery .....	26
6.2	Technical data, Environmental requirements .....	26
<b>7</b>	<b>Operation</b> .....	<b>27</b>
7.1	Control lights .....	27
7.2	Performing a scanning procedure .....	27
<b>8</b>	<b>Operational errors and faults</b> .....	<b>28</b>
<b>9</b>	<b>Cleaning</b> .....	<b>29</b>
<b>10</b>	<b>Returning old devices</b> .....	<b>29</b>
10.1	Disposing of the device .....	29
<b>11.</b>	<b>Addendum – Warning Signs</b> .....	<b>30</b>
	<b>EC conformity to type</b> .....	<b>31</b>

# 1 Indications for use

## 1.1 Indications for use

The Cercon® eye is a device for the **scanning of plaster models** as a manufacturing basis for bridges and caps in the **dental technology area**.

The device functions only with a Windows-based computer on which the Cercon® art software is installed and via which the Cercon® eye is controlled. The scan process occurs through laser scanning in conjunction with a special camera technology.

## 1.2 About this instruction for use

Observation of this operating manual is essential for proper and safe operation of the Cercon® eye. This operating manual contains important information for the safe, proper and efficient operation of the device. Observing this information helps to avoid further risks, repair costs and downtimes and helps to increase the reliability and operating lifetime of the Cercon® eye.

The operating manual must be read and fully understood by the Cercon® eye operator.

DeguDent GmbH assumes no liability for damages that arise due to improper operation of the Cercon® eye and/or disregarding the specifications of this operating manual.

## 1.3 Organisational elements used in this manual

Safety information, bodily injury, accident or damage warnings:

The device may only be opened by authorised DeguDent GmbH personnel!

### Step-by-step instructions:

1. Remove ...
2. Position ...

### Listed items:

1. Rotary plate ...
2. Laser outlet ...

### Identifier for operation and control elements:

**START** button

### In bold:

The green **control lamp (1)** goes out during the scan ...

## 1.4 Device type and year of manufacture

Cercon® eye, Year of manufacture: 2006

## 1.5 Manufacturer

DeguDent GmbH  
A Dentsply International Inc.  
Rodenbacher Chaussee 4  
D-63457 Hanau  
Phone +49 6181 59 5759

### Distributor

Dentsply International Inc.  
Dentsply Ceramco  
570 West College Avenue  
York, PA 17405-0872  
Phone (800) 877-0020

## 1.6 Copyrights and property rights

© October 2006, DeguDent GmbH

All rights to this operating manual are reserved, specifically the right of reproduction and distribution as well as the translation. In particular, this operating manual may not be reproduced or partially reproduced in any form (print, photo-copy, microfilm or by any other process) and/or stored in electronic systems without the prior approval of the DeguDent GmbH.

Violation is prohibited by law, and may result in civil and criminal penalties. DeguDent GmbH reserves all industrial property rights.

### 2.1 Contraindications

A different use or a use beyond the Cercon eye's scope does not comply with the regulations. The user/operator of the Cercon® eye is liable for any damages resulting thereof. This equally applies to unsanctioned changes to the device. Sanctioned use requires specific compliance with the information:

- on safety,
- on the operation,
- on maintenance, servicing and fault remedy measures that are described in this operating manual.

The device is only suitable for use in dental technology areas. For other areas of application or purposes, a prior written approval of DeguDent GmbH is required.

DeguDent GmbH devices are constructed and built using **state-of-the-art technology** and approved safety-technology regulations.

During use, **risks for personnel or other people or damages** to the Cercon® eye and other material assets can, may, result when the device:

- is operated by untrained or unauthorised personnel,
- is used in an unapproved way,
- is operated and serviced in an unapproved way.

Only allow **qualified and/or trained personnel** to perform the activities described in this manual. The legal age requirements for operating the device must also be observed!

Personnel in training may only operate the Cercon® eye **under constant supervision** of an experienced specialist or skilled employee!

### 2.2 Transport safety information

First visually check the **“Shockwatch”-Sticker (1)** on the carton. If its color has changed to red, the impact energy during transport was higher than permitted and your device could be damaged. Have your shipper confirm the triggering of the “Shockwatch” in writing.

Due to the device weight of approx. **31 kg/68.34 lb**, special precautions are required when transporting the Cercon® eye:

- **Open the packaging and carefully unpack the Cercon® eye.**
- **Use suitable load-handling equipment and load securing devices, if required.**
- **Only lift the device with at least two people!**
- **Check the device for external damage before setting up. A damaged device may not be set up and put into operation. In particular, the device may not be used if the plastic hood is damaged. Please contact us if your Cercon® eye is damaged.**

## 2 Contraindications

### 2.3 Operating safety information

To ensure a safe operation of the device, observe the following information:

- Before using Cercon® eye for the first time, a DeguDent accredited service technician will recheck the calibration to ensure accurate scan results. Please make sure that the device has been kept at the optimum working temperature (environmental temperature) of  $20\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$  ( $68\text{ }^{\circ}\text{F} \pm 2\text{ }^{\circ}\text{F}$ ) for at least 24 hours before the calibration is checked.
- After the calibration has been checked, Cercon® eye is not to be moved. A recheck of the calibration will become necessary if the device is moved in any way.
- Protect device from moisture and water vapour!
- If the device is not to be used for a longer period of time, disconnect the Cercon® eye from the mains supply!



## 3 Warnings

### 3.1 Warnings

The laser used within the digitalisation unit is a Class 2 laser with a maximum output of 2 mW. **Staring directly into the beam** is to be avoided with lasers of this power class; a momentary exposure (0.25 seconds), as would occur when accidentally looking into the beam is not considered dangerous. The use of optical corrective lenses (magnifying glass, binoculars) can increase the risk of eye damage. The following safety information related to the device is to be observed:

- Operating the Cercon® eye without detached housing is forbidden. This can lead to eye damage!
- Do not manipulate the turn-off device on the hood!
- Replace the warning information (position recognizable in the graphic) where necessary with the supplied warning information for your language.
- A closed hood does not prevent laser light from exiting the scanner. Do not look into the laser light even when the hood is closed!
- A defective hood must be replaced immediately for your own safety!

## 4.1 Precautions

The following information must be observed for maintenance and repair:

- **Maintenance and repair work may only be performed by technical personnel authorised by DeguDent GmbH!**
- **Only original replacement parts from DeguDent GmbH are to be used.**

The Cercon® eye is a device for **scanning plaster models** as a manufacturing basis for **bridges and caps** in dental technology areas. The device functions only with a Windows-based **computer (Cercon® art PC)** on which the **Cercon® art** software is installed and via which the Cercon® eye is controlled. The scan process occurs through laser projection in conjunction with two stereoscopically arranged **cameras**.

The models are mounted on a **modelling table (supplied)** and placed on the **rotary plate** located in the interior of the device. The precise **alignment** is supported via the Cercon® art software preview function. An **additional camera** supplies the image information for the adjustment. Lifting/closing the hood switches the **lighting** on/off automatically. The **scan process** is **started via the software after the hood has been closed**.

**Each tooth** is scanned **individually**. For bridges the scanning process (open hood, adjust position of abutment tooth or pontic, close hood, start scan-process) must be repeated for each unit until the complete bridge has been scanned.

The data or the **graphic models** are **displayed** within the Cercon® art software and can be **edited, saved** and **transferred** to other memory devices.

- **All further informations on working with the software can be found in the Cercon® art Instruction for use.**

## 4.2 Front side


The mechanical and electronic components of the Cercon® eye are protected by a solid metal **housing (1)**.

The **base unit (2)** consists of a aluminium plate to prevent vibrations during the scan.

The device is switched on via the **power switch (4)**.

The working area is covered by a pivoting, **transparent and darkened plastic hood (3)** which provides protection from the laser light and enables a restricted monitoring during the scan process. It can be opened by applying **light manual pressure** on the lower edge.

**Control lights (5)** are arranged above the working area.


## 4.3 Rear side


**The device may only be opened by authorised DeguDent GmbH personnel!**

The rear device opening is screwed together with a solid **covering plate (6)**; the following connections are located below.

- **Power connection (7.1), 100 – 240 V with integrated mains fuse (250 V, 1AT, high breaking capacity, 5 x 20 mm/5 x 0.78 in).**
- **USB 2 interface to the PC (7.2), for the communication with the PC/Cercon® art software.**

## 4 Precautions


### 4.4 Interior with rotary plate

A **rotary plate (8)** powered by a motor is located in the interior for positioning and aligning the modelling table. The interior  **matt black**  coating prevents light reflections that could cause false scan results.

In addition, there is an **opening (11)** in the upper area for the camera for alignment of the model (preview). The **lighting** occurs through this opening via a circle of diodes arranged around the camera.


The **laser openings (9)** as well as the openings for the stereoscopically arranged **cameras (10.1 and 10.2)** are located at the back of the chamber.

### 4.5 Modelling table with gauge marks

The **base plate (9.1)** on the **modelling table (9)** is used in order to bring the **plaster model (10)** into a firm scanning position. With the **clamping screw (9.2)** the model is held to this base plate, the adjustment is carried out via the **adjustment screw (9.3)** (hydraulic clamping device). The plaster teeth to be scanned should be aligned so that they lie at the approximate height of the gauge marks (9.4) and tooth removal direction is perpendicular to the modelling table.

**The position of the plaster model must not be changed during the entire scanning process but should rather be removed as a complete unit with the modelling table.**

The exact position of the modelling table is determined by the camera via the three **gauge marks (9.4)**. If the model is moved between the individual scan procedures (for multiple scans/bridges), the camera recognises the new position via the gauge marks.

The **gauge marks** must **never be covered** and must be **clearly visible** for the upper **camera (11)** during the entire scan. Dust, contamination or damages to the gauge marks can lead to poor multiple scans (e.g. poor bridge fitting).

### 5.1 Advers reactions

Not known

### 5.2 Setting up


1. Place the Cercon® eye on a **solid, flat and vibration free surface**. This must be suitable for absorbing the residual vibrations that occur during the scanning process. Otherwise **the working results** may be affected **adversely**.

2. **Strong light** can interfere with the cameras during the scanning process and **adversely affect the scan result**. The device should therefore not be operated in direct sunlight or near any other strong light source.

### 5.3 Electrical connection

3. Plug the **power cable (7.1)** into the device and connect it to a fused outlet (see connection load under **6.2** Technical data, Environmental requirements).

4. Connect the Cercon® eye and the computer using the **USB cable (7.2)**.

The initial setup and operation of the device may only be performed by authorised DeguDent GmbH personnel!

## 6 Scope of delivery

### 6.1 Included in the delivery

- Cercon® eye
- Model table (1)
- Scan holder (2)
- Scan template (3)
- Calibration ball (4)
- USB cable
- Power cable (2 versions; please request other version when needed)
- CD-ROM with the calibration data (please keep in safe custody)

### 6.2 Technical data, Environmental requirements

**Power supply:**

100 – 240 V AV, 50/60 Hz

**Permitted voltage fluctuations:**

not greater than  $\pm 10 \%$

**Max. power consumption:**

100 W

**Dimensions (width x height x depth):**

490 mm/19.30 in x 447 mm/17.6 in x 569 mm/22.4 in

**Weight:**

31 kg/68.34 lb

**Use:**

For indoor use only

**Min/Max ambient temperature:**

5 °C to 40 °C / 41 °F to 104 °F

In order to obtain valid scan-results, an ambient temperature range of 18 °C to 25 °C/64.4 °F to 77 °F, optimum 68 °F  $\pm$  2 °F, must be maintained.

Relative humidity 5–80 %

The Cercon® eye is put into operation by switching on the **power switch (4)** (see chapter 5.3 Electrical connection) and starting the associated PC. To get a good scan result switch on Cercon® eye 30 min. before starting the scan. The **modelling table with the attached plaster model** is positioned on the **rotary plate**. Good scan results can be reached with ochre-, sand- or other pastel colored gypsums.

After closing the **hood** and entering the scan command, the scan process starts automatically. For bridges, this procedure is repeated for the corresponding number of teeth to be scanned.

## 7.2 Performing a scanning procedure

During the scan, do not touch the Cercon® eye. The plastic hood must remain closed! For recording multiple teeth (bridges), the base plate (9.1) must not be adjusted. Only move the model holder on the rotary plate!

Do not switch off the device during the scan. This leads to errors in the geometric alignment of the teeth when conducting multiple scans!


All further informations on working with the software can be found in the Cercon® art Instruction for use.

## 7.1 Control lights

If the Cercon® eye is ready for operation, the **green control light (5a)** lights up.

During scanning, the **red control light (5b)** lights up. **The hood (4) must remain closed** during the scan.

If the **hood** is **opened** during the scan, the scanning process is terminated immediately. A new scanning procedure must be started via the software; the data of the aborted scan is lost (see also chapter 8 Operational errors and faults).

## 8 Operational errors and faults

Control panel/error	Account/measure
Device does not work with power switch turned on.	Power cable and/or USB cable not connected. Internal fuse has blown. <b>Fault repair may only be performed by authorised DeguDent GmbH personnel!</b>

The user is responsible for the proper return of the Cercon<sup>®</sup> eye if necessary. Please use the original packaging. Please keep the packaging for this eventuality. DeguDent GmbH assumes no liability for damages to the device due to improper packaging!

## 9 Cleaning

Clean the Cercon® eye **housing** and the **hood** with a **dry** or **lightly moistened cleaning cloth** as needed. Never use a wet cleaning cloth since **moisture** can **damage** electronic parts inside the device.

The **interior** of the device and the modelling table can be cleaned with a **soft brush** as needed. Pay attention that the matt coated interior surfaces is preserved. **Keep the hood closed** to avoid dirt getting into the scanning chamber. The gauge marks on the model holder must remain undamaged.

## 10 Returning old devices

### 10.1 Disposing of the device (for EU countries only)

The device is an electronic device according to the “**Act Governing the Sale, Return and Environmentally Sound Disposal of Electrical and Electronic Devices**”.


It was identified in accordance with the existing law and provided with this symbol.

The device is **not intended for private use**. It is manufactured and delivered for commercial use and is to be **disposed** by the **end user** according to the specifications of the Electrical and Electronic Equipment Act.

## 11 Addendum – Warning Signs

Replace **the warning information** (position recognizable in the graphic) where necessary with the supplied warning information for your language.

## EC conformity to type

DeguDent GmbH

P. O. Box 1364 · 63403 Hanau

GERMANY

www.degudent.com

We hereby declare the conformity of the manufactured item:

## Cercon<sup>®</sup> eye

Device for scanning of plaster models in the dental technology area.

With the fundamental requirements of the following EC directives:

1. EC Directive 98/37/EC (CE Mark) for Machinery
2. EC-Low Voltage Directive 73/23/EWG – Security Regulations for Electrical Fittings with a Rated Voltage between 50 V and 1000 V.
3. EC Directive 89/336/EWG – Electromagnetic Compatibility

The following harmonised norms were applied:

DIN EN ISO 12100-1: 2003; DIN EN ISO 12100-2: 2003;

EN 60204-1: 1997; EN 61010-1: 2001; EN 60825-1: 2001;

EN 61000-4-2; EN 61000-4-3; EN 61000-4-4; EN 61000-4-5; EN 61000-4-6;

EN 61000-4-11; EN 55011; EN 55011; EN 61000-3-2; EN 61000-3-3

Hanau, March 9, 2006

Dr. Udo Schusser

Research and development

Hanau, March 9, 2006

Torsten Schwafert

Head of dental technology business department

The text and images for this operating manual were prepared with great care. However, typographical or other incorrect specifications may occur. Please observe that DeguDent GmbH accepts no liability for typographical or specification errors.

<b>1</b>	<b>Indications for use</b> .....	<b>34</b>
1.1	Indications for use .....	34
1.2	About this instruction for use .....	34
1.3	Organisational elements used in this manual .....	34
1.4	Device type and year of manufacture .....	34
1.5	Manufacturer, Distributor .....	34
1.6	Copyrights and property rights .....	34
<b>2</b>	<b>Contraindications</b> .....	<b>35</b>
2.1	Contraindications .....	35
2.2	Transport safety information .....	35
2.3	Operating safety information .....	36
<b>3</b>	<b>Warnings</b> .....	<b>36</b>
3.1	Warnings .....	36
<b>4</b>	<b>Precautions</b> .....	<b>37</b>
4.1	Precautions .....	37
4.2	Front side .....	37
4.3	Rear side .....	37
4.4	Interior with rotary plate .....	38
4.5	Modelling table with gauge marks .....	38
<b>5</b>	<b>Advers reactions</b> .....	<b>39</b>
5.1	Advers reactions .....	39
5.2	Setting up .....	39
5.3	Electrical connection .....	39
<b>6</b>	<b>Scope of delivery</b> .....	<b>40</b>
6.1	Including in the delivery .....	40
6.2	Technical data, Environmental requirements .....	40
<b>7</b>	<b>Operation</b> .....	<b>41</b>
7.1	Control lights .....	41
7.2	Performing a scanning procedure .....	41
<b>8</b>	<b>Operational errors and faults</b> .....	<b>42</b>
<b>9</b>	<b>Cleaning</b> .....	<b>43</b>
<b>10</b>	<b>Returning old devices</b> .....	<b>43</b>
10.1	Disposing of the device (EU) .....	43
<b>11.</b>	<b>Addendum – Warning Signs</b> .....	<b>44</b>
	<b>EC conformity to type</b> .....	<b>45</b>

# 1 Indications for use

## 1.1 Indications for use

The Cercon® eye is a device for the **scanning of plaster models** as a manufacturing basis for bridges and caps in the **dental technology area**.

The device functions only with a Windows-based computer on which the Cercon® art software is installed and via which the Cercon® eye is controlled. The scan process occurs through laser scanning in conjunction with a special camera technology.

## 1.2 About this instruction for use


Observation of this operating manual is essential for proper and safe operation of the Cercon® eye. This operating manual contains important information for the safe, proper and efficient operation of the device. Observing this information helps to avoid further risks, repair costs and downtimes and helps to increase the reliability and operating lifetime of the Cercon® eye.

The operating manual must be read and fully understood by the Cercon® eye operator.

DeguDent GmbH assumes no liability for damages that arise due to improper operation of the Cercon® eye and/or disregarding the specifications of this operating manual.

## 1.3 Organisational elements used in this manual

Safety information, bodily injury, accident or damage warnings:

### Step-by-step instructions:

1. Remove ...
2. Position ...

### Listed items:

1. Rotary plate ...
2. Laser outlet ...

### Identifier for operation and control elements:

**START** button

### In bold:

The green **control lamp (1)** goes out during the scan ...

## 1.4 Device type and year of manufacture

Cercon® eye, Year of manufacture: 2006

## 1.5 Manufacturer

DeguDent GmbH  
A Dentsply International Inc.  
Rodenbacher Chaussee 4  
D-63457 Hanau  
Phone +49 6181 59 5759

### Distributor

Dentsply International Inc.  
Dentsply Ceramco  
570 West College Avenue  
York, PA 17405-0872  
Phone (800) 877-0020

## 1.6 Copyrights and property rights

© October 2006, DeguDent GmbH

All rights to this operating manual are reserved, specifically the right of reproduction and distribution as well as the translation. In particular, this operating manual may not be reproduced or partially reproduced in any form (print, photocopy, microfilm or by any other process) and/or stored in electronic systems without the prior approval of the DeguDent GmbH.

Violation is prohibited by law, and may result in civil and criminal penalties. DeguDent GmbH reserves all industrial property rights.

### 2.1 Contraindications

A different use or a use beyond the Cercon eye's scope does not comply with the regulations. The user/operator of the Cercon® eye is liable for any damages resulting thereof. This equally applies to unsanctioned changes to the device. Sanctioned use requires specific compliance with the information:

- on safety,
- on the operation,
- on maintenance, servicing and fault remedy measures that are described in this operating manual.

The device is only suitable for use in dental technology areas. For other areas of application or purposes, a prior written approval of DeguDent GmbH is required.

DeguDent GmbH devices are constructed and built using **state-of-the-art technology** and approved safety-technology regulations.

During use, **risks for personnel or other people or damages** to the Cercon® eye and other material assets can, may, result when the device:

- is operated by untrained or unauthorised personnel,
- is used in an unapproved way,
- is operated and serviced in an unapproved way.

Only allow **qualified and/or trained personnel** to perform the activities described in this manual. The legal age requirements for operating the device must also be observed!

Personnel in training may only operate the Cercon® eye **under constant supervision** of an experienced specialist or skilled employee!

### 2.2 Transport safety information

First visually check the **“Shockwatch”-Sticker (1)** on the carton. If its color has changed to red, the impact energy during transport was higher than permitted and your device could be damaged. Have your shipper confirm the triggering of the “Shockwatch” in writing.

Due to the device weight of approx. **31 kg/68.34 lb**, special precautions are required when transporting the Cercon® eye:

- **Open the packaging and carefully unpack the Cercon® eye.**
- **Use suitable load-handling equipment and load securing devices, if required.**
- **Only lift the device with at least two people!**
- **Check the device for external damage before setting up. A damaged device may not be set up and put into operation. In particular, the device may not be used if the plastic hood is damaged. Please contact us if your Cercon® eye is damaged.**

## 2 Contraindications

### 2.3 Operating safety information

To ensure a safe operation of the device, observe the following information:

- Before using Cercon® eye for the first time, a DeguDent accredited service technician will recheck the calibration to ensure accurate scan results. Please make sure that the device has been kept at the optimum working temperature (environmental temperature) of  $20\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$  ( $68\text{ }^{\circ}\text{F} \pm 2\text{ }^{\circ}\text{F}$ ) for at least 24 hours before the calibration is checked.
- After the calibration has been checked, Cercon® eye is not to be moved. A recheck of the calibration will become necessary if the device is moved in any way.
- Protect device from moisture and water vapour!
- If the device is not to be used for a longer period of time, disconnect the Cercon® eye from the mains supply!


## 3 Warnings

### 3.1 Warnings

The laser used within the digitalisation unit is a class IIIa laser (CFDA, CFR 1040.10) with a maximum output of 2 mW. Staring directly into the beam or view directly with optical instruments (magnifying glass, binoculars) is to be avoided with lasers of this power class; The following safety information related to the device is to be observed:

- Operating the Cercon® eye without detached housing is forbidden. This can lead to eye damage!
- Do not manipulate the turn-off device on the hood!
- Replace the warning information (position recognizable in the graphic) where necessary with the supplied warning information for your language.
- A closed hood does not prevent laser light from exiting the scanner. Do not look into the laser light even when the hood is closed!
- A defective hood must be replaced immediately for your own safety!

## 4.1 Precautions

The following information must be observed for maintenance and repair:

- **Maintenance and repair work may only be performed by technical personnel authorised by DeguDent GmbH!**
- **Only original replacement parts from DeguDent GmbH are to be used.**

The Cercon® eye is a device for **scanning plaster models** as a manufacturing basis for **bridges and caps** in dental technology areas. The device functions only with a Windows-based **computer (Cercon® art PC)** on which the **Cercon® art** software is installed and via which the Cercon® eye is controlled. The scan process occurs through laser projection in conjunction with two stereoscopically arranged **cameras**.

The models are mounted on a **modelling table (supplied)** and placed on the **rotary plate** located in the interior of the device. The precise **alignment** is supported via the Cercon® art software preview function. An **additional camera** supplies the image information for the adjustment. Lifting/closing the hood switches the **lighting** on/off automatically. The **scan process** is **started via the software after the hood has been closed**.

**Each tooth** is scanned **individually**. For bridges the scanning process (open hood, adjust position of abutment tooth or pontic, close hood, start scan-process) must be repeated for each unit until the complete bridge has been scanned.

The data or the **graphic models** are **displayed** within the Cercon® art software and can be **edited, saved** and **transferred** to other memory devices.

- **All further informations on working with the software can be found in the Cercon® art Instruction for use.**

## 4.2 Front side


The mechanical and electronic components of the Cercon® eye are protected by a solid metal **housing (1)**.

The **base unit (2)** consists of a aluminium plate to prevent vibrations during the scan.

The device is switched on via the **power switch (4)**.

The working area is covered by a pivoting, **transparent and darkened plastic hood (3)** which provides protection from the laser light and enables a restricted monitoring during the scan process. It can be opened by applying **light manual pressure** on the lower edge.

**Control lights (5)** are arranged above the working area.


## 4.3 Rear side


**The device may only be opened by authorised DeguDent GmbH personnel!**

The rear device opening is screwed together with a solid **covering plate (6)**; the following connections are located below.

- **Power connection (7.1), 100 – 240 V with integrated mains fuse (250 V, 1AT, high breaking capacity, 5 x 20 mm/5 x 0.78 in).**
- **USB 2 interface to the PC (7.2), for the communication with the PC/Cercon® art software.**

## 4 Precautions


### 4.4 Interior with rotary plate

A **rotary plate (8)** powered by a motor is located in the interior for positioning and aligning the modelling table. The interior  **matt black**  coating prevents light reflections that could cause false scan results.

In addition, there is an **opening (11)** in the upper area for the camera for alignment of the model (preview). The **lighting** occurs through this opening via a circle of diodes arranged around the camera.


The **laser openings (9)** as well as the openings for the stereoscopically arranged **cameras (10.1 and 10.2)** are located at the back of the chamber.

### 4.5 Modelling table with gauge marks

The **base plate (9.1)** on the **modelling table (9)** is used in order to bring the **plaster model (10)** into a firm scanning position. With the **clamping screw (9.2)** the model is held to this base plate, the adjustment is carried out via the **adjustment screw (9.3)** (hydraulic clamping device). The plaster teeth to be scanned should be aligned so that they lie at the approximate height of the gauge marks (9.4) and tooth removal direction is perpendicular to the modelling table.

**The position of the plaster model must not be changed during the entire scanning process but should rather be removed as a complete unit with the modelling table.**

The exact position of the modelling table is determined by the camera via the three **gauge marks (9.4)**. If the model is moved between the individual scan procedures (for multiple scans/bridges), the camera recognises the new position via the gauge marks.

The **gauge marks** must **never be covered** and must be **clearly visible** for the upper **camera (11)** during the entire scan. Dust, contamination or damages to the gauge marks can lead to poor multiple scans (e.g. poor bridge fitting).

### 5.1 Advers reactions

Not known

### 5.2 Setting up


1. Place the Cercon® eye on a **solid, flat and vibration free surface**. This must be suitable for absorbing the residual vibrations that occur during the scanning process. Otherwise **the working results** may be affected **adversely**.

2. **Strong light** can interfere with the cameras during the scanning process and **adversely affect the scan result**. The device should therefore not be operated in direct sunlight or near any other strong light source.

### 5.3 Electrical connection

3. Plug the **power cable (7.1)** into the device and connect it to a fused outlet (see connection load under **6.2** Technical data, Environmental requirements).

4. Connect the Cercon® eye and the computer using the **USB cable (7.2)**.

The initial setup and operation of the device may only be performed by authorised DeguDent GmbH personnel!

## 6 Scope of delivery

### 6.1 Included in the delivery

- Cercon® eye
- Model table (1)
- Scan holder (2)
- Scan template (3)
- Calibration ball (4)
- USB cable
- Power cable (2 versions; please request other version when needed)
- CD-ROM with the calibration data (please keep in safe custody)

### 6.2 Technical data, Environmental requirements

#### Power supply:

100 – 240 V AV, 50/60 Hz

#### Permitted voltage fluctuations:

not greater than  $\pm 10 \%$

#### Max. power consumption:

100 W

#### Dimensions (width x height x depth):

490 mm/19.30 in x 447 mm/17.6 in x 569 mm/22.4 in

#### Weight:

31 kg/68.34 lb

#### Use:

For indoor use only

#### Min/Max ambient temperature:

5 °C to 40 °C / 41 °F to 104 °F

In order to obtain valid scan-results, an ambient temperature range of 18 °C to 25 °C/64.4 °F to 77 °F, optimum 68 °F  $\pm 2$  °F, must be maintained.

Relative humidity 5–80 %

The Cercon® eye is put into operation by switching on the **power switch (4)** (see chapter 5.3 Electrical connection) and starting the associated PC. To get a good scan result switch on Cercon® eye 30 min. before starting the scan. The **modelling table with the attached plaster model** is positioned on the **rotary plate**. Good scan results can be reached with ochre-, sand- or other pastel colored gypsums.

After closing the **hood** and entering the scan command, the scan process starts automatically. For bridges, this procedure is repeated for the corresponding number of teeth to be scanned.

## 7.2 Performing a scanning procedure

During the scan, do not touch the Cercon® eye. The plastic hood must remain closed! For recording multiple teeth (bridges), the base plate (9.1) must not be adjusted. Only move the model holder on the rotary plate!

Do not switch off the device during the scan. This leads to errors in the geometric alignment of the teeth when conducting multiple scans!


All further informations on working with the software can be found in the Cercon® art Instruction for use.

## 7.1 Control lights

If the Cercon® eye is ready for operation, the **green control light (5a)** lights up.

During scanning, the **red control light (5b)** lights up. **The hood (4) must remain closed** during the scan.

If the **hood** is **opened** during the scan, the scanning process is terminated immediately. A new scanning procedure must be started via the software; the data of the aborted scan is lost (see also chapter 8 Operational errors and faults).

## 8 Operational errors and faults

Control panel/error	Account/measure
Device does not work with power switch turned on.	Power cable and/or USB cable not connected. Internal fuse has blown. <b>Fault repair may only be performed by authorised DeguDent GmbH personnel!</b>

The user is responsible for the proper return of the Cercon® eye if necessary. Please use the original packaging. Please keep the packaging for this eventuality. DeguDent GmbH assumes no liability for damages to the device due to improper packaging!

## 9 Cleaning

Clean the Cercon® eye **housing** and the **hood** with a **dry** or **lightly moistened cleaning cloth** as needed. Never use a wet cleaning cloth since **moisture** can **damage** electronic parts inside the device.

The **interior** of the device and the modelling table can be cleaned with a **soft brush** as needed. Pay attention that the matt coated interior surfaces is preserved. **Keep the hood closed** to avoid dirt getting into the scanning chamber. The gauge marks on the model holder must remain undamaged.

## 10 Returning old devices

### 10.1 Disposing of the device (for EU countries only)

The device is an electronic device according to the “**Act Governing the Sale, Return and Environmentally Sound Disposal of Electrical and Electronic Devices**”.


It was identified in accordance with the existing law and provided with this symbol.

The device is **not intended for private use**. It is manufactured and delivered for commercial use and is to be **disposed** by the **end user** according to the specifications of the Electrical and Electronic Equipment Act.

## 11 Addendum – Warning Signs

Replace **the warning information** (position recognizable in the graphic) where necessary with the supplied warning information for your language.

## EC conformity to type

DeguDent GmbH

P. O. Box 1364 · 63403 Hanau

GERMANY

www.degudent.com

We hereby declare the conformity of the manufactured item:

## Cercon<sup>®</sup> eye

Device for scanning of plaster models in the dental technology area.

With the fundamental requirements of the following EC directives:

1. EC Directive 98/37/EC (CE Mark) for Machinery
2. EC-Low Voltage Directive 73/23/EWG – Security Regulations for Electrical Fittings with a Rated Voltage between 50 V and 1000 V.
3. EC Directive 89/336/EWG – Electromagnetic Compatibility

The following harmonised norms were applied:

DIN EN ISO 12100-1: 2003; DIN EN ISO 12100-2: 2003;

EN 60204-1: 1997; EN 61010-1: 2001; EN 60825-1: 2001;

EN 61000-4-2; EN 61000-4-3; EN 61000-4-4; EN 61000-4-5; EN 61000-4-6;


EN 61000-4-11; EN 55011; EN 55011; EN 61000-3-2; EN 61000-3-3

Hanau, March 9, 2006

Dr. Udo Schusser

Research and development

Hanau, March 9, 2006

Torsten Schwafert

Head of dental technology business department

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Het apparaat mag alleen door geautoriseerd vakpersoneel van DeguDent GmbH worden geopend!