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Instructions for Use Solaris Flash-Liquid

(GB)

Solaris Flash-Liquid is a highly acidic and cyanidic gold bath to be used as a preparatory gold-plating bath for friction renewal and for gold-plating cast-metal dental objects. Solaris Flash-Liquid in the Solaris unit lets you precipitate an extremely thin and homogenous film of gold on metal surfaces previously activated with Solaris Pre-Liquid in a matter of only a few minutes. Solaris Flash-Liquid ensures the best possible bond between the base to be gold-plated and the subsequent electroplated gold layer of 5 to 30 µm in thickness applied using Solaris gold baths.

Amount of liquid 500 ml
Movement in bath None

Gold content 1 g/500 ml
Colour of liquid Clear pink

Working temperature Room temperature
Colour of precipitate Shiny golden

Contraindications: Do not use in known hypersensitivity to gold or cobalt.

Warnings: May cause cancer by inhalation. Harmful by inhalation, in contact with skin and if swallowed. Contact with acids liberates very toxic gas. Irritating to eyes and skin. Hydrocyanic acid may cause all levels of intoxication.

Protective measures: Avoid aerosol and steam formation (use neither steam cleaners or pressurized air for cleaning and/or drying). Use your personal protective gear (protective goggles, gloves). For other important information on protective measures, consult the Material Data Safety Sheet (MSDS).

Adverse effects: Reactions that may occur are allergic reactions against metals contained in the liquid or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the liquid have been reported in isolated cases. Sensitization through skin contact is a possibility.

Interactions: Contact with acids (including carbonic acid/carbon dioxide) will release hydrocyanic acid that is flammable and may form explosive gas/air mixtures. Sensitization through skin contact is a possibility. Avoid occlusal and proximal contacts of different metals.

Storage: Avoid contamination and direct sunlight. After use, keep container tightly closed and in a well-ventilated place, protected from light.

Recommend storage temperature: 5-35 °C (41-95 °F).

Use before date: See label

Processing notes:

Dosage: Use 500 ml of liquid in a large beaker. If required, add more Solaris Flash-Liquid until all objects to be gold-plated are covered.

Caution: Solaris Flash-Liquid must never be mixed with Solaris gold baths. Avoid contamination, and thoroughly rinse the beakers under running water after use.

For additional processing instructions using the Solaris unit, consult the "Solaris - Instructions for Use for Friction Renewal and Gold-Plating". Pre-plating is performed in the Solaris unit at room temperature and without stirring. Solaris Flash-Liquid is suitable for multiple use. It can be used between 20 and 30 times or until the colour transition during preparatory pre-plating becomes less pronounced. This presupposes that the liquid is stored in the closed original container after use and protected from light.

Disposal considerations:

Product: Dispose of in accordance with all applicable laws and regulations.

Recommendation: Residue or unusable solutions should be turned over to a recognized waste management company. Package and label all waste as if it were the pure substance. Do not remove the labels on shipping containers until they are finally disposed of.

Uncleaned containers:

Do not reuse empty containers. Dispose of in accordance with all applicable laws and regulations. When handling "empty" containers containing product residue, follow the instructions given on the containers' labels.

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