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Degutron[®] eco

Instructions for Use
Instructions d'utilisation
Istruzioni per l'uso
Manual de instrucciones

D	Gebrauchsanweisung	4-15
GB	Instructions for Use	16-27
F	Instructions d'utilisation.....	28-39
I	Istruzioni per l'uso	40-51
E	Manual de instrucciones.....	52-63

D

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The text and illustrations of these Instructions for Use have been compiled with the utmost diligence. Nevertheless, the presence of typographic errors or incorrect data cannot be excluded. Please note that DeguDent GmbH will not be responsible for such errors.

- ➔ **Content..... 17**

- 1. Introduction..... 18–20**
 - 1.1. General description of the Degutron® eco..... 18
 - 1.2. Machine components..... 18
 - 1.3. Scope of delivery..... 18
 - 1.4. Technical Data..... 18
 - 1.5. General information..... 19

- 2. Putting into operation..... 21**
 - 2.1. Set-up directions..... 21
 - 2.2. Mains supply..... 21
 - 2.3. Cooling water..... 21
 - 2.4. Vacuum..... 21
 - 2.5. Air pressure supply..... 21

- 3. Operation..... 22–24**
 - 3.1. Operation elements..... 22
 - 3.2. Melting and casting..... 23

- 4. Service..... 25**
 - 4.1. Trouble shooting..... 25
 - 4.2. Maintenance..... 25

- 5. EC Declaration of Conformity..... 26**

- 6. Error and Warning Codes..... 27**

1 Introduction

1.1. General description of the Degutron® eco

With this casting machine, the melting charge is heated up by induction. Following the liquid material is cast into the mould.

The special advantage of inductive heating is that the melting temperature is reached rapidly, because the heat is directly generated in the metal. During the melting process, the metal is thoroughly mixed by means of a magnetic field, ensuring a homogenous mixing.

1.2. Machine components

The machine consists of separate units built into one casing:

- a) Inside the casing, the following components are installed:
- microprocessor controlled induction generator
 - medium frequency transformer
 - oscillation circuit capacities
- b) The front panel consists of:
- operating panel for control of melting cycles

- c) The casting unit consists of:

- inductor
- crucible
- insulation parts

- d) The backside covers:

- RS 232 Port for service
- Connectors for water, vacuum, air pressure

1.3. Scope of delivery

One each unless otherwise noted:

- Vacuum pressure casting unit
- Instructions for Use
- Short Operating Instructions
- Silica crucible shield
- Ceramic crucible
- Ceramic crucible with graphite insert
- Transparent display foil (gold)
- Protective plug for thermal element socket
- 2 pressure hoses, 6 × 1 mm, 2 m long
- 2 pressure hoses, 8 × 1.25 mm, 2 m long
- Compressor connector for 6-mm tube
- 2 pressure connectors QS-1/4-6
- Pressure connectors QS-1/4-8
- Quartz fleece for crucible shield
- Protective material for chamber floor
- 4 muffle spacers (1x, 3x, 6x, 9x)
- Vacuum pump

1.4. Technical Data

Crucible volume	15 cm ³ (ceramic crucible)
	10 cm ³ (with graphite inlay)
Crucible temperature	max. 2000 °C
Melting performance	3.5 kW
Mains supply	230 V / 16 A, 50 or 60 Hz
Cooling water supply	2.5–5 bar / min. 1.5 l/min
Cooling water recoil	Pressureless
Cooling water entry temperature	10–20 °C
Air pressure supply	4–6 bar
Vacuum	min. 8 m ³ /h, 0.1 mbar abs.
Room temperature	10–35 °C
Relative atmospheric humidity	20–80%
Weight	ca. 27 kg
Dimensions (width x depth x height)	400 mm x 400 mm x 450 mm

1.5. General information

Safety information

In order to ensure a constant, ideal performance of the machine and to ensure safe working conditions, the user is to observe the following safety measures:

- The complete electrical wiring is to be performed only by qualified and specially trained personnel.
- Check the machine and the supply devices regularly on possible damages.
- Do not reach into the lock mechanism while closing.
- The top of the pressure chamber has to be positioned horizontally.
- When opening cabinet/casings or when removing parts, certain parts under electric power can cause danger. If opening up the machine is necessary (before maintenance, change of machine settings, repair or exchange of parts) the machine must be disconnected from mains supply. If working on the opened machine is inevitable, only by DeguDent authorized, qualified trained personnel aware of the danger caused hereby and aware of the relevant regulations may be instructed to do this work.
- Capacitors in the machine can still be charged even when the machine's mains supply is switched off.
- When it seems that the machine can no longer be used safely it has to be taken out of use and secured against further unintentional use. The following incidents indicate that safe working is no longer possible:
 - The machine is visibly damaged.
 - The machine does not function.
 - Uncommonly heavy wear from transportation.
- The safety valves (pressure relieve valves) placed inside the machine casing may not be removed, closed or altered in any way.
- As melted metal is processed with this machine (= temperatures up to 2000 °C), the utmost care and attention has to be applied when working with the machine. The following necessary safety clothing is principally to be worn:
 - fireproof clothes
 - fireproof closed shoes
 - fireproof gloves
 - safety goggles
- Special caution is essential when working with graphite crucibles and graphite dies because these are only **visibly** hot at temperatures above 500 °C!
- In commercial enterprises the regulations for the prevention of accidents of the relevant authority for electrical machinery are to be followed.
- People with pacemakers must not be near the machine while it is running.
- A trained responsible personnel is to supervise work with this machine.

1 Introduction

Use as determined/liability

This casting machine is designed and manufactured to comply with the latest technical knowledge and according to approved safety regulations. When not putting into operation correctly or not using the machine as determined, however, danger and damages may occur. Therefore we suggest reading carefully and completely this manual before putting the machine into operation and to follow the instructions given therein:

- This machine is only to be used for melting, and casting of common dental alloys. Melting unknown alloys can be hazardous by toxic substances. Its released exhaust air shall not be set free. The risk analysis shall be evaluated in advance.
- Do not melt substances which are inflammable or can form an explosive gas during melting.
- For melting already inserted alloys appropriate safety cloths are required such as rubber gloves, masks, etc.
- This machine may only be connected as mentioned in this manual. The mains supply as well as in- and outgoing pressures are to be provided as stated on the machine label.
- This machine was designed for use in closed spaces (indoors) and may only be used for the above-mentioned purpose.
- This unit is designed for a maximum of six castings per hour.
- DeguDent does not take over responsibility for any damage, caused by non original DeguDent consumables and spare parts, unless the damage occurred independently.
- The machine may not be changed or modified. Technical changes may be accomplished only with the written permission of DeguDent GmbH.
- Disregarding or wrong interpretation the data in this operating instructions DeguDent is not responsible for damage, unless the damage occurred independently.
- This instruction manual is conform to the latest technical condition of the machine when printed. Technical changes and fittings are subject to change.

On delivery, immediately check if the machine is complete and if there are any transportation damages.

2.1. Set-up directions

The machine should be placed on clean and dry ground. The ground should be even, solid and level. The rear side of the machine must remain freely accessible to ensure sufficient circulation of cooling air. The cooling air temperature must not exceed 35 °C and must not be contaminated. Only after all supply and connecting systems have been connected, the machine may be put into operation.

2.2. Mains supply

Only trained personnel may connect the machine. Please pay attention to the nominal voltage and frequency stated on the machine label. The 1-phase current may differ +/- 10% (maximum) from the nominal value. The mains supply must at least be furnished with 16 A delay-action fuses (at 1* 230 V AC). The machine is equipped with a 16 A-CeeKon-connector when leaving INDUTHERM. All 3 lines (L1, N, PE) must be connected correctly. The wall socket must be freely accessible, so that the equipment can be shut down in emergency by pulling the power supply.

2.3. Cooling water

Use two hoses with an outer diameter of 6 mm for the cooling water. The water entry pressure must be at least 2.5 bar and must not be more than 5 bar. The water drainage must be without pressure. The water entry temperature must be above 10 °C and must not exceed 20 °C. The concentration of lime in the water must not be more than 60 mg CaO/l. The cooling water must not be contaminated.

The cooling water supply must be open as long as it is still a hot crucible inside the machine.

2.4. Vacuum

The vacuum pump should be connected with the supplied tube with an outer diameter of 8 mm. Keep the vacuum tube as short as possible to avoid pressure loss.


2.5. Air pressure supply

Use a hose with an outer diameter of 8 mm for air pressure. The entry pressure must be 4–6 bar. Keep the air tube as short as possible to avoid pressure loss. Avoid during the casting process additional air pressure load to ensure the appropriate pressure. Otherwise the casting quality can be negatively influenced.

3 Operation

3.1. Operation elements

The upper display shows the “Power/Temperature/Time” target values.

The lower display ‘Pressure’ shows the current pressure values and switches to timer mode (seconds).

Upper Display “Power/Temperature/Time”

“---”	Initial state “ready”
“P000 to P100”	Heating power in % (10% steps), or In tilted position after casting display shows a timer in seconds.
“0.50 to 2.00”	Casting pressure, displayed after pressing “+” and “-” simultaneously. The factory-provided pressure is 2.00 bars.

“Exxx” Error code

Buttons “+” und “-”

-	To set the heating power in %
-	To set the casting pressure (press “+” and “-” simultaneously, set value with “+” or “-” between 0.50–2.00 bars) for parameter setting
-	If “+” and “-” is not further pressed, the display changes back to “ready”.

Button „Start“

1x press Heating on and evacuation of casting chamber

Repeated press Change between full vacuum and reduced vacuum (please see respective recommendations of the alloy producer). This function is activatable only after reaching the target vacuum level.

Button “Stop”

1x press Running program will be aborted, generator off, the vacuum in the casting chamber will be released to ambient pressure.

Lower Display “Pressure”

Shows the actual pressure (absolute) inside the casting chamber “-0.10 to -1.00” vacuum

“xxxx” If a vacuum of < -0.98 bars is reached, the display switches to a timer counting in seconds. The maximum possible vacuum is achieved after approx. 30 seconds.

“0.00” Ambient pressure

“0.00 to 2.00” Overpressure

Service Connection (RS232)

Connection for service/diagnostics at the backside of the machine

3.2. Melting and Casting

Before starting check crucible and crucible insulation for dirt residues/possible damages.

Insert a small piece of the white quartz fleece into the crucible shield.

The crucible should touch slightly the quartz fleece, but the collar of the crucible should align with the top of the crucible shield.

Crucible shield with crucible in the induction coil.

3 Operation

General casting process:

- Open cooling water supply.
- Open air pressure supply.
- Switch on vacuum pump.
- Switch on main switch.
- After successful self test the display „Power“ shows “---”. (ready)
- Before every first casting the machine performs a short pressure test. The casting chamber is set under a slight overpressure before the vacuum starts. If this test fails (e.g. the casting chamber lid is not closed properly, an error code is displayed, e.g. E081)

To reduce a temperature loss of the mould, this test can be preceded as follows:

- Insert empty crucible with graphite inlay with the quartz tube into the induction coil.
 - Close and lock casting chamber.
 - Press “Start” 1x to run the test. The test is successfully completed if the chamber is getting evacuated and the heating starts.
 - End test by pressing “Stop” 1x. Insert the desired crucible. **Caution, the graphite crucible is hot!**
-
- Insert crucible with material into the quartz glass tube inside the inductor.
 - Place the respective mould holder into the casting chamber. The holder can be fixed to keep moulds of 50–80 mm length in 5 mm steps.
 - With “+” or “-” the desired heating power from 0–100% can be adjusted. This is also possible during melting.
 - Place mould in the mould holder.
 - Close and lock casting chamber.
 - Press “Start” 1x, unit automatically switches to vacuum inside casting chamber.
 - If the achieved vacuum is higher than -0.98 bar, the display switches to a timer.
 - By repeated pushing the ‘Start’ button you can choose between full and reduced vacuum.

- If casting temperature is reached (see data sheet of manufacturer of the alloy) press the button on the handle to release and tilt the unit uninterrupted clockwise till the catch. Release the button on the handle to lock the unit in the tilted position. During the tilting the unit automatically switches from vacuum to overpressure inside the casting chamber. The generator stops. On the display “Time” a timer counts the seconds after casting.
- After a sufficient waiting time (min. 1 min.), unlock the handle and tilt the unit back to horizontal position. The vacuum inside the crucible chamber is released to ambient pressure.
- Open the casting chamber and remove the mould and the mould holder.

Important note for the use of the ceramic crucible without graphite:

The induction system needs for a good coupling into the metal a certain mass/volume of metal.

If you use metal ingots, place the ingots horizontal, parallel to the crucible bottom. Also short cylinders should be placed this way for a better heating performance. Leave enough space to the crucible wall, because the heated metal will expand and can break the crucible.

For the specific use of graphite or ceramic crucibles please follow the instructions of the respective alloy producer.

4.1. Trouble shooting

Only trained personnel should open the machine.

The machine cannot be turned on:

- when there is no mains supply

The heating cannot be turned on:

- when there is no cooling water supply
- any error code „Exxx“, see appendix 5.2.

4.2. Maintenance

Caution: Make sure to disconnect unit from mains before any service or maintenance!

Check crucibles for damage before use and replace as appropriate.

DeguDent GmbH
Rodenbacher Chaussee 4
63457 Hanau
Germany

28 May 2009

We hereby declare that the product:

Degutron® eco

Induction-heated vacuum casting unit for melting precious and non-precious metal alloys

is in conformity with the basic provisions of the following Directives:

1. Directive 2002/95/EC on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)
2. Directive 2002/96/EC on waste electrical and electronic equipment (WEEE)
3. Directive 2006/95/EC, low-voltage directive
4. Directive 2004/108/EC, EMC directive

The following standards were applied:

DIN EN 61010-1:2001

DIN EN 62076:2006

DIN EN 55011:2007 and A2:2007

DIN EN 61000-6-2:2005

DIN EN 61000-6-3:2007

Hanau, 28 May 2009
Dr Udo Schusser
Director, Research and Development

Hanau, 2 June 2009
Dr Bernhard Kraus
Director, Quality Management

Error/ Warning	Description	Cause	Action
E012	Water pressure low	<ul style="list-style-type: none"> • Tubes too long • Pressure too low • Dirt 	<ul style="list-style-type: none"> • Keep the connections as short as possible • Water pressure should be 2.5–5.0 bar with min. flow of 1.5 l/min • Clean circuit and exchange tubes
E021	Over temperature power unit	<ul style="list-style-type: none"> • Unit can not cool down • Too many consecutive casts • Water cooling insufficient 	<ul style="list-style-type: none"> • Do not block the space at the back of the machine • Not more than 6 casts/hour • See E012
E030	IZK Peak current too large	<ul style="list-style-type: none"> • Not enough material in the crucible • Old crucible • No crucible inserted 	<ul style="list-style-type: none"> • Only heat with inserted filled crucible
E040	Generator switched on/off too often	<ul style="list-style-type: none"> • Max. 5 switches/min 	<ul style="list-style-type: none"> • Switch off the casting unit and avoid too many switching events
E050	Over temperature generator	<ul style="list-style-type: none"> • Unit can not cool down • Too many consecutive casts • Water cooling insufficient 	<ul style="list-style-type: none"> • Do not block the space at the back of the machine • Not more than 6 casts/hour • See E012
E051	Power on time too long	<ul style="list-style-type: none"> • Too many casts per hour 	<ul style="list-style-type: none"> • Not more than 6 casts/hour
E055	No vacuum	<ul style="list-style-type: none"> • Lit not closed and locked • No vacuum connected or dysfunctional • Lit lock not tight 	<ul style="list-style-type: none"> • Close and lock lit • Control lit ring for damages or dirt • Check vacuum pump and connections • Tighten lit lock
E064	Power on time too long	<ul style="list-style-type: none"> • Too many casts per hour 	<ul style="list-style-type: none"> • Not more than 6 casts/hour
E065	Generator power too low	<ul style="list-style-type: none"> • More than 15% deviation 	<ul style="list-style-type: none"> • Check electric power supply
E066	Generator power too high	<ul style="list-style-type: none"> • More than 15% deviation 	<ul style="list-style-type: none"> • Check electric power supply
E081	Casting chamber pressure out of tolerance	<ul style="list-style-type: none"> • Lit not locked • No sufficient pressure supply • No vacuum supply • Lit lock not tight 	<ul style="list-style-type: none"> • Close and lock lit • Control lit ring for damages or dirt • Check pressure supply (4–6 bar) • Check vacuum pump and connections • Tighten lit lock

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