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Motorcast compact

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

DeguDent

A Dentsply International Company

Motorcast compact

D	Gebrauchsanweisung	5
GB	Instructions for use	11
F	Mode d'emploi	17
I	Istruzioni per l'uso	23
E	Instrucciones para el uso	29

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Motorcast compact

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Motorcast compact

Scope of delivery

The scope of delivery of the Motorcast compact includes:

- Basic unit, Motorcast compact
- Centrifugal arm for torch melting
- Instructions for use
- Melting crucible

Accessories available

- Universal melting crucible **REF 2351 0054**
- Transferable control board **REF 2351 1502**
- Automatic Gas-saving MC **REF 2343 0001**

Technical Data

Voltage for operation:	See identification plate
Max. power loading:	250 W
Required utility protection:	See identification plate
Operating Temperature:	4 °C bis 40 °C
Overall dimensions:	475 mm high, 505 mm wide, 550 mm depth
Weight, incl.centrifugal arm:	ca. 410 N (42 kg)

Safety precautions

Read these instructions carefully before installation and start-up. Always follow the advice given in these instructions.

The Motorcast compact is to be run strictly with original replacement parts and accessories. Only in this way can the rated operational loadings be achieved and the required safety assured in the workplace.

Dangerous voltages are present within the unit. Maintenance work and repairs are to be performed only by authorized DeguDent service technicians. Protect the equipment from dampness. This will avoid short circuiting and keep the unit free from damage.

Check the correctness of the seating and functional safety of the centrifugal arm before each operation.

If the unit is not going to be used for a long time, unplug the power cord.

Instruction for the user

The Motorcast compact operates on the centrifugal casting principle.

It is designed for the DeguDent casting ring system. The use of gypsum and phosphate-based investments is possible, as well as the use of ring less casting systems.

The Motorcast compact is equipped with an electronic, centrifugal, RPM controller, which permits casting speed to be varied in a range between 350 to 450 rpm. Thus, the unit can be adapted to the investment employed, alloy types, and characteristics of the dental casting.

The RPM setting can be adjusted on the control display of the panel board.


Figure 1: Control Display

Initial Operation - Placement of the unit

Place the Motorcast compact on a stable, level stand, elevated to working height.

As with any centrifugal casting equipment, the apparatus will shake slightly with every run. For this reason, insure that equipment is on a stable, level stand so that it will not slide off.

Installation and removal of the centrifugal arm


Figure 2: Securing the Centrifugal Arm

To remove the centrifugal arm, remove the centrally located fastening screw and lift the arm upward.

To install, place the centrifugal arm, as indicated in Fig. 2, on the drive shaft, then replace the central fastening screw and turn this down tightly.

Be careful that the two carrier rods of the drive shaft are securely seated in the drillings provided in the centrifugal arm.

Adjustment of the counterweight

Before each casting, the counterweight must be adjusted for the size of the ring employed. Only in this way can a smooth run be achieved.

Poorly set counterweights can create excess vibrations and thus cause the unit to slip.

Fixed positioning on the rotating centrifugal arm must be in accord with the DeguDent casting ring system (1x, 3x, 6x, 9x, 12x). To adjust the counterweight, hold down the lock button and slide the weight to the desired position. Take care that the counterweight is seated in the selected position, i.e. the lock button protrudes from the counterweight.

A counterweight which is not locked in position can be displaced, leading to strong vibrations, which can cause the unit to slip.

In the secured position, a more refined adjustment of the counterweight is possible by raising and turning the outermost adjustment ring (see Fig. 2). Make sure that all parts are once again locked into the preset positions.

The settings +1 and +2 are available for heavy rings requiring a larger mass of alloy and settings -1 and -2 for lighter rings and smaller quantities of alloys. These may be read on the central axis of the centrifugal arm.

Casting with the centrifugal arm for torch melting

Select an RPM commensurate to the cast object. The rotation speed may be set between 350 and 450 RPM. When a phosphate investment is used we recommend a speed of 450 RPM. When using gypsum investment however, because of lower compressive strengths, a lower RPM, for example 350 RPM, should be used. To set the rotational speed, use the buttons on the control board. The preset speed will be indicated by means of an illuminated diode on the button bar.

Rotate the centrifugal arm in its enclosure in such a manner that the rising, hot gases do not accumulate inside the housing.

The ring, the crucible slides, as well as other parts are subject to extreme heating when open flame is used for melting. Moreover, the housing can heat up depending on the position of the centrifugal arm and the burners, and on how long the flame burns.

Parts that can withstand high temperatures during operation should never be touched unless your hands are protected. Burns can result.

During the first start-up, odors may issue from the centrifugal arm and the housing insulation. Discoloration of the housing insulation is also possible.

Place the casting crucible in its position and correct, if necessary, the height of the ring support to match that of its casting mold. Do this by means of the height adjustment on the centrifugal arm.

The slide carriage for the casting crucible is provided with an automatic locking mechanism, which prevents the crucible from touching the ring during the casting and thus damaging the latter.

To adjust the slides, lift the adjustment yoke, which also releases the locking mechanism, and then push the slides into the desired position. The adjustment lever is designed so that it can easily be handled with the tongs for the ring.

Following the melting of the alloy, the casting device is made ready for casting by pushing the "Pour" button. When it is ready to pour, a signal is given by a yellow pilot light. The actual casting will begin when the housing is closed.

In order to avoid damage to the centrifugal arm, casting may only take place when the casting ring is in place and the counterweight is properly adjusted. Make sure that there are no foreign objects in the casting space.

Lay the preheated ring on its support. Lift the adjustment lever off the crucible slide with the mold tongs and push the casting crucible all the way to the mold. Be sure that the casting crucible lines up exactly with the casting trough of the ring. Possible differences in height between the casting crucible and the ring can be corrected by turning the height adjustment on the centrifugal arm.

Before the release of the casting process, shove the casting crucible against the mold. Otherwise both the casting crucible and the mold can be damaged when the casting process begins.

The control sequence for the initiation of the casting process can be changed, i.e. by first closing the housing and then starting the process by pushing the “CAST” button. In this case, casting begins immediately after the CAST button is pushed.

Conclude the casting process by pushing the CAST button again. The yellow control light will be extinguished and the centrifugal arm will stop. The housing safety lock will be released immediately after the arm comes to rest.

Altering position of control board

With the auxiliary “control conversion set”, the control board can be moved from its position at the rear pedestal to the left side of the unit.

This change enables you to have the installation of the Motorcast compact under a low shelf, extending outward.

Such changes should only be made by DeguDent service personnel.

Installation of the Automatic Gas-saving MC

The new Automatic Gas-saving MC is designed especially for use with the Motorcast compact.

It is outstanding because it has a high degree of operating safety, is always ignition-ready, and can be adjusted for a predetermined flame characteristic. In addition, it has the advantages of easy manipulation and economy in gas use. It is available as a separate option and can be placed on the Motorcast compact with very few attachments.


Fig. 3: Automatic Gas-saving device

Care and Maintenance

Care and maintenance comprise the activities which keep the Motorcast compact in an impeccable operating condition.

Examine the unit regularly for mechanical damage. Ascertain that the seating of the centrifugal arm as well as the counterweighting is in order, as well as that of the housing itself. In this way, you will avoid accidents and downtime on the machine.

Thoroughly remove any accumulations of contaminants on moving parts. By doing so, you are assuring easy operation and satisfactory functioning of the equipment.

Before cleaning the unit, turn off the casting apparatus switch and also unplug the power cord.

When cleaning the machine, be careful to prevent any moisture from penetrating it. This will remove the danger of short circuits which can lead to malfunction of the machine or fire. Dangerous voltages are present in the inner wiring of the equipment. Maintenance and repair work should therefore be carried out only by authorized DeguDent service technicians.

Making changes to the unit yourself lessens the operational safety of the casting machine and increases the danger of personal injury and material damage.

Trouble Shooting - Tripped circuit-breakers/ blown fuses

Upon start-up or during running of the equipment, the external building circuit-breaker or fuse protection shuts off current.

The on-rush of current upon start-up of the motor is many times greater than the normal operating current. Check whether or not too many users are on the line which serves the casting machine. Preheating ovens, for example, are great consumers of energy.

Upon start-up or during running of the equipment the internal overload protection trips.

The protective device of the equipment is found beside the power switch, and current can be reestablished by pressing the reset button. Repeated drop out indicates a defect in the casting machine which must be corrected by an authorized service technician from DeguDent.

Declaration of Conformity

DeguDent GmbH

P.O. Box 13 64
D-63403 Hanau

We certify herewith the conformity of the product:

Motorcast compact/ centrifugal casting apparatus

with the basic requirements of the following EC-Directives:

- 1. EC-Directive for machines 98/392/EWG including change 1, 91/368/EWG, Change 2, (93/44/EWG), Change 3, 93/68/EWG**
- 2. EC-Directive for electrical machinery for use within specified ranges of voltage (low-voltage directive) 73/23/EWG with a Change 1, 93/68/EWG**
- 3. EC-Directive for electro-magnetic compatibility 89/336/EWG with a Change 1, 93/68/EWG**

The following associated standards were also applied:


**EN 292-1: 1991, EN 292-2: 1991, EN 60204-1: 1992, EN 61010-1: 1993,
EN 50081-1: 1991, EN 50082-1: 1991**

The following international technical standards were observed:


IEC 801-1: 1984, IEC 801-2: 1984

The following national technical standards were observed:

UVV VBG 7z: 4/81

Torsten Schwafert
Director business area
dental laboratory products

Dr. Alexander Völcker
Research & Development

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