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Multicast compact

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

Inhaltsverzeichnis / Table of contents / Sommaire / Indice / Indice

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Instructions for use

Multicast compact

Electrical power-driven melting and casting unit for dental alloys.

1 Working with the Multicast compact

1.1 Scope of delivery

The scope of delivery of Multicast compact includes:

- basic unit Multicast compact
- centrifuge arm for electric melting
- high-performance melting furnace HOMC
- operating instructions
- mirror
- graphite crucible

1.2 Accessories available

Centrifuge arm for flame melting	REF 2351 1401
Universal melting crucible	REF 2351 0054
High-performance melting furnace HOMC	REF 2351 1482
Graphite crucible	REF 2351 0068
Ceramic crucible	REF 2351 0069
Conversion kit control element	REF 2351 1402
Multicast Printer	REF 2351 1403
Spare paper roll for Multicast Printer	REF 2351 1404
Automatic gas-saving unit MC	REF 2343 0001

1.3 Technical data

Power supply:	Voltage see typeplate
Maximum power consumption:	Heating power 1.25 kW Drive power 0.25 kW
Mains fuse required:	See typeplate
Ambient temperatures:	4 °C to 40 °C
Dimensions (h x w x d):	475 mm x 505 mm x 550 mm
Weight incl. centrifuge arm:	approx. 510 N (52 kg)

1.4 Safety information

Carefully read these operating instructions prior to installation and putting into operation. The information contained in these operating instructions must be observed.

The Multicast compact must only be operated with original spare parts and accessories listed. Only in this way the performance data given can be achieved and the required operational safety is ensured.

Dangerous voltages are applied in the inside of the unit. Therefore maintenance and repair work may only be carried out by service technicians authorized by DeguDent.

Protect the unit against moisture in order to avoid short circuits and damage to the unit.

Before the unit is activated, the correct fixation and proper function of the centrifuge arm must always be checked.

Unplug the unit if it is not used over longer periods.

1.5 Information for the user

The Multicast compact works according to the principle of centrifugal casting. It is designed as a modular system. The basic unit is equipped with an electric melting device. With the accessory „centrifuge arm for flame melting" which is additionally available the Multicast compact can be easily converted to melting with the open flame.

The Multicast compact is suitable to receive the DeguDent casting ring system. Processing of plaster- and phosphate-bound investment compounds as well as the use of flasks are possible.

The Multicast compact has a powerful microprocessor for operator prompting, motor control and control of the electric melting furnace. User entries are possible with the key pad of the control element. Messages and help information are displayed in the integrated LC-display ①.


Figure 1: key pad

- ① = LC-display
- ② = numerical block
- ③ = readiness for casting

The casting speeds can be varied within a range of 350 rpm to 450 rpm. Due to that, the unit can be easily adjusted to the individual investment compounds, alloy types and dental technical objects.

2 Putting the unit into operation

2.1 Installation of the unit

Place the Multicast compact on a flat, stable surface.

As it is common for centrifugal casting machines, a unit slightly jerks when it is started. Therefore the unit must be placed on a flat, stable surface to avoid that the unit will get out of place.

Exchanging the centrifuge arm

With the accessory „centrifuge arm for flame melting" that is available separately the Multicast compact can be easily converted to melting with the open flame.

For this purpose remove the central fixation screw with the hexagon key that is supplied together with the „centrifuge arm for flame melting".

Remove the centrifuge arm by pulling it upward and place the new centrifuge arm on the drive shaft as shown in figure 2.

It must be ensured that both driving pins of the drive shaft engage in the respective borings of the centrifuge arm.

Insert the central fixation screw again and tighten it with the help of the enclosed hexagon key.


Figure 2: Fixation of the centrifuge arm

2.3 Adjustment of the counterweight

Prior to each casting process the counterweight must be adjusted to the casting ring size that is used. Only in this way smooth running of the machine is ensured.

A counterweight that is not correctly adjusted may causes strong vibrations and leads to the fact that the unit gets out of its original position.

The lock-in positions on the centrifuge arm are matched with the DeguDent casting ring system (1x, 3x, 6x, 9x). For readjustment of the counterweight press in the locking knob and slide the weight to the desired position with the locking knob pressed in. Please ensure that the counterweight engages in the desired position, i.e. the locking knob projects again from the counterweight.

A counterweight that is not locked may shift its position in an uncontrolled way, cause strong vibrations and lead to the fact that the unit gets out of its original position.

In the locked position fine adjustment of the counterweight can be carried out by lifting and turning the outer adjusting ring (see figure 2). It must be ensured that the adjusting ring engages in the preselected position i.e. the locking knob projects from the counterweight.

For heavier casting rings or large alloy quantities the positions +1 and +2 must be selected; for lighter casting rings or smaller alloy quantities select the positions -1 and -2. Reading is performed at the center axle of the centrifuge arm.

2.4 Casting with the centrifuge arm for electric melting

Turn on the unit with the power switch on the right side of the unit. The display ① will briefly show the serial number of your unit. Over the display in the control element you will now be requested to enter the furnace temperature.

Enter the desired furnace temperature with the keys of the numerical block ②. The furnace temperature can be adjusted between 700 °C and 1450 °C.

If, e.g., you wish to have a furnace temperature of 900 °C, enter the figures 0900. After pressing the last key, the entry is automatically completed. You may also confirm the value in the display by pressing the „E" key.

If you have made an error you can change the current entry at any time by pressing the „C" key, or return to the previous input field after completing the entry.

After preselecting the furnace temperature a further fusion time in seconds can be entered. With the activation of the casting readiness the further fusion time is automatically started with the „Gießen (cast)" key ③ after the alloy has melted. It can be selected in a range from 0 to 255 seconds.

When the entry of the further fusion time has been completed, the casting speed is requested. Preselect a speed that is suitable for the casting object. The speed can be adjusted between 350 rpm and 450 rpm. In the case of phosphate-bound investment compounds we recommend a speed of 450 rpm. Due to the lower compressive strength of plaster-bound investment compounds, a lower speed, e.g. 350 rpm should be used.

When the temperature, further fusion time and casting speed have been entered, the microprocessor will check the correct position of the centrifuge arm. If the centrifuge arm is not in the contact position required for heating up, the display and an acoustic signal will request you to turn the centrifuge arm to the contact position. The arm can be turned freely.

Immediately after the centrifuge arm has engaged in its position, the furnace is heated up to the preadjusted temperature. The current furnace temperature is shown in the display.

Furnace, casting ring support as well as adjacent parts will heat up strongly during the heating process of the furnace. Do not touch these parts after the unit has been started because of the risk of getting burned.

When heating up the unit for the first time, the furnace, centrifuge arm and housing insulation may produce bad smells.

Shortly before the preselected furnace temperature is reached, a short acoustic signal can be heard which indicates the readiness of the furnace.

Insert your graphite or ceramic crucible into the furnace now. The furnace flap remains closed during this process. Before you fill in the alloy, wait until the crucible is completely warmed up (graphite crucible 5 min, ceramic crucible approx. 10 min).

Adjust the height of the casting ring support so that it is suitable for your casting ring. Adjusting is performed with the adjusting knob on the side of the centrifuge arm.

After the alloy has melted, the casting unit is prepared for casting by pressing the „Gießen (cast)” key ③. The readiness for casting is shown in the display. The casting process can be started by closing the housing.

In order to avoid damage to the furnace and the centrifuge arm, casting must only be carried out if the casting ring has been inserted and the counterweight is adjusted correctly. In addition to that, make sure that there are no objects in the casting chamber when the casting process is started.

The „Gießen (cast)” key ③ activates the readiness for casting and simultaneously starts the clock for the further fusion time. The expired time is indicated in the display. An acoustic signal will inform the user about the end of the preselected further fusion time. The acoustic signal can be deactivated by pressing the „E” key.

Place the preheated casting ring on the casting ring support, open the flap of the melting furnace completely and slide the furnace with the adjusting lever of the furnace to the casting ring so that the crucible base projects into the funnel of the casting ring. Any differences in height can be adjusted by turning the height adjustment at the centrifuge arm. The adjustment lever of the furnace has been designed so that it can be easily operated with casting ring tongs. Afterwards, move the crucible e.g. with the crucible tongs, entirely against the casting ring.

Before starting the casting process, furnace and crucible must be moved entirely against the casting ring. Otherwise the furnace, casting ring and crucible might be damaged when the casting process is started.

After the end of the further fusion time, close the housing in order to start the casting process.

The further fusion time can be reset by pressing the „NEW” key. Casting can be started by closing the housing independent of the expiration of the further fusion time.

The operating sequence to start the casting process can also be changed, i.e. first close the housing and then start casting by pressing the „Gießen” (cast) key ③. In this case casting is started immediately after pressing the „Gießen” (cast) key ③.

End the casting process by pressing the „Gießen” (cast) key ③ again. The yellow indicator light goes out and the centrifuge arm is stopped. The lock of the housing opens automatically immediately after the centrifuge arm stops.

In case the casting process is not switched off manually, the unit will automatically stop after 5 minutes.

2.5 Use of the program selection keys

Settings of temperature, further fusion time and speed that are more frequently used can be programmed freely (see chapter 3.3) and can be fetched with program selection keys.

Select a preset program by pressing the program key „PROG” followed by the program number (1-9). The display will briefly show the stored program name. The furnace temperature, further fusion time and speed are automatically preselected.

The display shows the programmed speed which you can still change with the numerical block so that it is suitable for your casting object. Confirm the program selected by pressing the „E” key.

2.6 Casting with the centrifuge arm for flame melting

With the accessory „centrifuge arm for flame melting” that is available separately the Multicast compact can be easily converted to melting with the open flame.

If required, the centrifuge arm must be converted as described in chapter 2.2.

Switch on the unit and enter the temperature value „0” to adjust the unit to the „flame melting operation mode”. During this operation mode the temperature is not requested again until the unit is switched off. Conversion to „furnace melting operation mode” can be performed after exchanging the centrifuge arm and switching on the unit again by entering a temperature higher than „0”.

Enter the required further fusion time and casting speed for your alloy and casting object (see chapter 2.4).

Turn the centrifuge arm so that rising, hot gases are not accumulated in the housing.

Casting ring support, crucible carriage as well as adjacent parts strongly heat up during melting with the open flame. Additionally, the housing may heat up depending on the position of the centrifuge arm, the adjustment of the burner as well as the burning time of the flame.

Do not touch parts that heat up during operation without suitable protection because of the risk of getting burned.

When heating up the unit for the first time, the centrifuge arm and housing insulation may produce bad smells. Discoloration of the housing insulation is possible.

Due to the high starting speed the special melting crucible (slotted type) must not be used since a back draft can occur in front of the discharge opening and therefore splashing of alloys can result. If required, readjust the height of the casting ring support with the height adjustment at the centrifuge arm.

The carriage for the melting crucible features an automatic locking device so that the melting crucible does not press against the casting ring and damage it.

For moving the carriage, lift the adjustment bar to release the locking device and move the carriage to the desired position. The adjustment lever is designed so that it can be easily operated with casting ring tongs.

In order to avoid damage to the centrifuge arm, casting must only be carried out if the casting ring has been inserted and the counterweight is adjusted correctly. In addition to that, make sure that there are no objects in the casting chamber when the casting process is started.

By pressing the „Gießen“ (cast) key ③ the clock for the further fusion time is started automatically. The expired time is indicated in the display. An acoustic signal will inform the user about the end of the preselected further fusion time. The acoustic signal can be deactivated by pressing the “E” key.

Place the preheated casting ring on the casting ring support. Lift the adjustment lever of the crucible carriage with the casting ring tongs and move the melting crucible entirely against the casting ring. Ensure that the melting crucible projects exactly into the funnel of the casting ring. Any differences in height can be adjusted by turning the height adjustment at the centrifuge arm.

Before starting the casting process, the melting crucible must be moved entirely against the casting ring. Otherwise the melting crucible and the casting ring might be damaged when the casting process is started.

After the end of the further fusion time, close the housing in order to start the casting process.

The further fusion time can be reset by pressing the „NEW“ key. Casting can be started by closing the housing independent of the expiration of the further fusion time.

The operating sequence to start the casting process can also be changed, i.e. first close the housing and then start casting by pressing the „Gießen“ (cast) key ③. In this case casting is started immediately after pressing the „Gießen“ (cast) key ③. The further fusion time will not run down.

End the casting process by pressing the „Gießen“ (cast) key ③ again. The yellow indicator light goes out and the centrifuge arm is stopped. The lock of the housing opens automatically immediately after the centrifuge arm stops.

In case the casting process is not switched off manually, the unit will automatically stop after 5 minutes.

2.7 Activating the clock

In order to print out date and time in the protocol of a connected documentation printer, the clock of the Multicast compact must be activated. It is activated by inserting a battery into the control element. The battery ensures that the clock is still working after the unit has been switched off.

For this purpose switch off your Multicast compact. Loosen the two fixation screws ④ of the control element and carefully remove the control element by pulling it upward. Insert a 1.5 V round cell (type AA) in the battery case of the control element. Observe the polarity indicated in the battery case.

Reinsert the control element and fasten it with the two fixation screws ④. Make sure that the connecting cable ⑤ is not wedged in when inserting the control element.

Set date and time as described in chapter changing date and time. After the activation of a documentation printer (see chapter 2.9) the display briefly shows the date and time when the unit is switched on.

2.8 Moving the control element

With the „conversion kit for the control element“ that is available separately the control element can be moved from the control pillar at the back to the left side of the unit.

The conversion allows you to use the Multicast compact even under a low rack system.

Switch off and unplug your Multicast compact before moving the control element.

1.) Loosen the two fixation screws ④ of the control element as shown in figure 3 and carefully remove the control element from the control pillar by pulling it upward. Pull off the plug of the connecting cable ⑤ at the control element and put the cable into the control pillar.


Figure 3


Figure 5

2.) Close the control pillar with a cover that is contained in the conversion kit ⑩. Make sure that the connecting cable ⑤ is not wedged in when fastening the cover. Screw the two screws ④ into the respective borings of the cover ⑩.


Figure 4


Figure 6

3.) Remove the fixation screws ⑥ of the cover plate at the bottom of the unit and replace the cover plate with the control element socket or the Multicast Printer. Reinsert the two fixation screws ⑥ and tighten them. Connect the plug of the connecting cable ⑦ to the control element and insert the control element into the control element holding device. Screw the two fixation screws ⑧ in the respective borings of the control element.

4.) Connect the cable ⑨ of the control element holding device with the ⑨ interface of the Multicast compact.

2.9 Use of the Multicast Printer


Figure 7: Printer ports

With the Multicast Printer that is available as an accessory you can document the most important casting parameters for quality assurance if the centrifuge arm for electric melting is used.

The casting protocols that are printed out contain the parameters of alloy, alloy quantity, batch, furnace temperature and casting speed. Additionally, the laboratory name, date and time are automatically entered into the casting protocol.

The Multicast Printer is connected with the Multicast compact via a special interface ⑨ and screwed firmly to the housing. Please observe the separate instructions for assembly. It is prepared to receive the control element (figure 5). It is also possible to connect commercially available, ASCII-compatible printers with a serial port to the Multicast compact via the integrated serial interface.

Prior to connecting the printer, switch off the Multicast compact and connect the printer to the casting unit. First switch on the printer and then the Multicast compact and select the external casting protocol as described in chapter 3.7.

Connection of an external printer to the Multicast compact unit

It is possible to connect an external printer to the Multicast compact unit. A printer with serial interface, e.g. HP Deskjet™, is required. The interface of the printer must be set - in accordance with the description in the printer manual - to the following transmission mode:

Serial, 9600 baud, 8 bit, not parity bit („none“), 1 stop bit, handshake via DTR, no Xon/ Xoff.

A cable with the following configuration must be used as a transmission cable:


Pin 9-pole socket (casting unit)	Pin 25-pole plug (printer)	Signal (casting unit)	Signal (printer)	Direction	Meaning
3	3	TxD	RxD	to printer	Printer data to printer
5	7	Gnd	Gnd	-	signal mass
8	20	CTS	DTR	from printer	printer feedback „ready to receive”

Putting the printer into operation

- 1.) Connect the printer with the 9-pole socket of the Multicast compact (the 25-pole socket must only be used to connect the original printer or the control element!)
- 2.) Switch on the printer for „online” mode
- 3.) Switch on the casting unit without pressing any key. Only in the case of initial recognition the message „Drucker erkannt (printer recognized)” will be displayed briefly. If this message is not displayed, check the cable once more and make sure that the printer is „online”.
- 4.) Switch off the casting unit and switch on again with the „PROG” key pressed down.
- 5.) Change the casting protocol to “extern (external)” in the „Grundeinstellungen vornehmen (perform basic settings)” mode. The corresponding menu will only appear, if the printer has been recognized correctly (see chapter 3.7). Exit „Grundeinstellungen vornehmen (perform basic settings)”.

Now, with the printer being activated („online”), the message „Bitte warten, Ausdruck (please wait, print-out)” should be displayed briefly after each casting process and a printout should be performed.

Troubleshooting diagram after recognition of printer

2.10 Use of the automatic gas-saving unit MC

The automatic gas-saving unit MC has especially been developed for the use with the Multicast compact.

It features high operational safety, tremendous readiness for ignition and constantly preselectable flame adjustment, economic consumption as well as comfortable handling. It is available as an additional accessory and can be easily mounted to the Multicast compact.

Note!

Gas installations must only be performed by authorized specialists (gas fitter).

3 Basic settings

3.1 General information

In order to switch over to the „Grundeinstellungen vornehmen (perform basic settings)” mode first switch off the Multicast compact. Press the „PROG” program key and hold it down while switching on the Multicast compact.

You are now in the „Grundeinstellungen vornehmen (perform basic settings)” mode. The display shows the language. The following basic settings can be carried out:

Figure 8: automatic gas-saving unit MC

- Changing the language
- Configuration of program keys 1 - 9
- Activating and deactivating the acoustic signal
- Changing date and time
- Changing the laboratory name
- Activating or deactivating the printer
- Adjustment of the melting furnace.

Press the „E” key several times in order to display the sequence of basic settings.

You can exit the „Grundeinstellungen vornehmen (perform basic settings)” mode with the „C” key. You will automatically get back to the standard operation mode.

3.2 Changing the language

Select „Sprache (language)” with the „E” key in the display as described in chapter 3.1. Switch between the languages with “NEW” until the desired language is shown in the display. It is possible to choose between German, English, French, Italian and Spanish. Confirm the language selected by pressing the „E” key or exit the program with the „C” key without making any changes.

3.3 Configuration of the program selection keys 1-9

Select the basic setting of the program to be configured with the „E” key as described in chapter 3.1. If, e.g., you would like to change the setting of program key 4, press the enter key („E”) until „P4” is shown at the bottom of the display. Press the „NEW” key to change the configuration of the program key.

The cursor is now directly in the input field for the program name. The letter at the cursor position can be changed now with the keys 2 and 8. The keys 4 and 6 will move the cursor within the input field.

Complete the entry of the program name with the „E” key. The cursor will automatically return to the input field for the temperature (T) of the respective program.

Accept the preadjusted temperature value by pressing the enter key or enter a new value with the numerical block of the key pad.

The respective speed (D) and further fusion time (W) for the program key can be entered in the same way.

Changes within an input field can be cancelled with the „C” key. The „NEW” key will move the cursor to the beginning of the input field.

Complete the entry of the program with the „E” key. The data are stored and the next setting is displayed.

3.4 Activating and deactivating the acoustic signal

Select the setting for the acoustic signal with the enter key „E” as described in chapter 3.1. Activate or deactivate the acoustic signal with „NEW”. When the acoustic signal is deactivated, no acoustic messages will be issued. Warning signals cannot be suppressed.

Confirm by pressing the „E” key or exit the program part with the „C” key without making any changes.

3.5 Changing date and time

Select the setting for date and time with the enter key „E” as described in chapter 3.1. Press „NEW” to change date and time.

Use the numerical keys to enter the date in the form of DD.MM.YY, e.g. 040504 for May 4, 2004. After entering the last figure the cursor will automatically return to the input field for the time.

The time is also entered with the numerical keys in the form of HH.MM, e.g. 1245 for 12:45.

Changes within an input field can be cancelled with the „C” key. The “NEW” key will move the cursor to the beginning of the input field.

Complete the entry of date and time with the enter key „E”. Date and time are changed and the next setting is displayed.

3.6 Changing the laboratory name

The laboratory name is only displayed in the casting protocol of a connected printer. It is only useful to enter the laboratory name if a documentation printer is used.

In order to enter the laboratory name, select the setting for the laboratory name with the enter key „E” as described in chapter 3.1. Press „NEW” to change the laboratory name.

The letter at the cursor position can be changed now with the keys 2 and 8. The keys 4 and 6 will move the cursor within the input field.

Complete the entry of the laboratory name with the enter key „E”. Changes within the input field can be cancelled with the “C” key.

3.7 Activating and deactivating the casting protocol

This program part is only displayed if a connected documentation printer has been recognized when the unit is switched on.

Select the setting for the casting protocol with the enter key „E” as described in chapter 3.1. With the „NEW” key you shift between „extern (external)”, „intern (internal)” and „aus (off)”. Select „extern” to send the casting protocol to a conventional, ASCII-compatible printer via the integrated serial interface. With the „intern” setting you direct the data output to the documentation printer which is available as an accessory. By selecting „aus (off)” the printout of the casting protocol is avoided.

Confirm your selection with the „E” key or exit the program part with the „C” key without making any changes.

3.8 Adjustment of the melting furnace

Select the setting for the furnace adjustment with the enter key „E” as described in chapter 3.1. Press „NEW” in order to change the set value.

Prior to delivery each melting furnace is subject to a test program in which a characteristic number is determined which considers the individual temperature distribution in the furnace. The characteristic number is enclosed to each furnace and should be set in order to ensure optimum temperature control and a long service life of the melting furnace.

The characteristic numbers for fine adjustment can be set within the range between 1 and 9. If you do not know the characteristic number of your furnace, please set the value 5.

Confirm your selection with the „E” key or exit the program part with the „C” key without making any changes.

4 Cleaning and maintenance

Cleaning and maintenance include all types of work which guarantee the proper working state of the Multicast compact.

Within regular intervals check the unit on mechanical damage and correct position of the centrifuge arm, the counterweight as well as the protection housing.

Immediately remove all contaminations from movable construction parts. This will ensure smooth running and proper function.

Prior to cleaning the unit must be switched off and unplugged.

During cleaning the unit make sure that no moisture reaches into the unit. This will avoid short circuits which may result in a breakdown or a fire.

In the inside of the unit dangerous voltages are applied. Maintenance work and repair must therefore only be performed by service technicians authorized by DeguDent.

Manipulations of the unit by the user reduce the operational safety and include the risk of damage to persons and materials.

5 Elimination of errors

5.1 Fuse breaks the circuit

During activation or start of the unit the external domestic fuse breaks the circuit.

The starting current exceeds the operating current by far. Therefore it must be checked whether too many devices are operated in the electric circuit of the casting unit. Preheating furnaces e.g. consume a lot of energy.

During activation or start of the unit the internal fuse breaks the circuit.

The unit fuse is located next to the power switch and is reset by pressing in the button. Breaking the circuit indicates a defect of the unit which must be eliminated by a service technician authorized by DeguDent.

5.2 Temperature display

During heating up of the furnace the temperature is not displayed correctly or the value varies strongly.

Remove the centrifuge arm as described in chapter 2.2. Check the contacts at the unit and the centrifuge arm. If required, clean the contacts with fine abrasive paper (grain size 400-800).

DeguDent GmbH

P.O. Box 1364
D-63403 Hanau

This is to certify the conformity of the product:

Multicast compact/centrifugal casting unit as of unit number 42001395

with the basic requirements of the following EC guidelines:

- 1. EC guideline for machines 89/392/EEC with 1. Modification 91/368/EEC, 2. Modification (93/44/EEC) and 3. Modification 93/68/EEC**
- 2. EC guideline for electrical equipment to be used within certain voltage limits (Guideline for Low Voltage) 73/23/EEC with 1. Modification 93/68/EEC**
- 3. EC guideline for electromagnetic compatibility 89/336/EEC with 1. Modification 93/68/EEC**

The following harmonized standards were applied:

EN 292-1: 1991, EN 292-2: 1991, EN 60204-1: 1992, EN 61010-1: 1993, EN 50081-1: 1991, EN 50082-1: 1991

The following international standards and specifications were applied:

IEC 801-1: 1984, IEC 801-2: 1984

The following national technical standard and specification was applied:

UVV VBG 7z: 4/81

The procedure was adhered to in accordance with Appendix V of the guideline for machines.

Dr. Josef Rothaut
Director (Manager) Research & Development

Erich Steiner
Director (Manager) operative unit

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