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# Cercon<sup>®</sup> ceram Kiss ht-Dentine

Ergänzung zur Gebrauchsanweisung  
Directions for Use – Supplement  
Complément au mode d'emploi  
Integrazione delle istruzioni per l'uso  
Ampliación de las instrucciones de uso

# Cercon® ceram Kiss ht-Dentine

<b>D</b>	Ergänzung zur Gebrauchsanweisung	4
<b>GB</b>	Directions for Use – Supplement	7
<b>F</b>	Complément au mode d'emploi	10
<b>I</b>	Integrazione delle istruzioni per l'uso	13
<b>E</b>	Ampliación de las instrucciones de uso	16

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## General note

Cercon® ceram Kiss ht-dentines supplement the Cercon® ceram Kiss range of veneering ceramics by providing the option of veneering the highly translucent Cercon® ht framework material using a simple dentine/incisal build-up in V shades A1 to D4. The general Directions for Use for Cercon® ceram Kiss apply.

## Purpose

- Cercon® ceram Kiss ht-dentines were developed specifically for veneering crowns and bridges made of Cercon® ht light and medium.
- For dental use only.

For information on contraindications, warning messages for medical devices, safety notes and technical specifications, please consult the general Directions for Use for Cercon® ceram Kiss.

## Instructions

A reproducible tooth shade for the veneered Cercon® ht restoration can be obtained by combining the appropriate framework shade (Cercon® ht light, Cercon® ht medium) with the corresponding dentine/incisal build-up (Table 1).

## Notes

The shade and appearance of highly translucent framework materials are influenced by the shade of the underlying prepared tooth. In addition, non-vital and/or severely discoloured hard tissue or coloured cementing agents may further affect the aesthetic result.

Moreover, highly variable cross-sections of the framework will also have an impact on shade. These effects can be compensated for by applying suitable ceramic stains.

**Table 1: Cercon® ceram Kiss ht-dentines shade-matching chart for V shades**

V shade	Cercon® ht framework	Dentine	Incisal
A1	Light	A1	1
A2	Medium	A2	2
A3	Medium	A3	3
A3,5	Medium	A3,5	3
A4	Medium	A4	5
B1	Light	B1	1
B2	Light	B2	1
B3	Medium	B3	4
B4	Medium	B4	6
C1	Light	C1	1
C2	Light	C2	5
C3	Medium	C3	5
C4	Medium	C4	6
D2	Medium	D2	2
D3	Medium	D3	4
D4	Medium	D4	4

- *Cercon® ceram Kiss is not compatible with any other veneering ceramics.*
- *Cercon® ceram Kiss must not be used together with Cercon® ceram love due to differences in material properties and firing temperatures.*

**Table 2: General firing recommendations – Cercon® ceram Kiss**

Firing	Pre-heating °C	Drying time min	Heating rate °C/min	Final temp °C	Holding time min	Vacuum hPa	Long-term cooling min
Dentine 1	450	5:00	55	830	1:30	50	–
Dentine 2	450	5:00	55	820	1:30	50	–
Glaze firing	450	3:00	55	800	1:00	–	6:00
Correction (Final-Kiss)	450	5:00	55	680	1:00	50	6:00
Final shoulder (F-SM)	450	5:00	55	680	1:00	50	6:00

Slow cooling is mandatory after the last firing; this includes correction firings of restorations after try-in. The firing temperature must be adapted to the number of units fired in the same cycle. Five to ten units require an increase by 5°C to 10°C; more than ten units require an increase by 10°C to 20°C. The values listed here are intended for orientation only and should be regarded only as guidelines. Your firing results may differ. All firing results depend on the performance of the furnace used, which in turn depends on the make, model and age of the furnace. Therefore, the guideline values will have to be adapted individually for each firing. We recommend running a test firing cycle to evaluate the performance of the furnace used. We have compiled and checked all values and other data with great care. However, we cannot be liable for your results under any circumstances.

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DeguDent GmbH