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Cercon[®] PressMaster

Instructions for use

Mode d'emploi

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DeguDent

A Dentsply International Company

Cercon[®] PressMaster

D	Gebrauchsanweisung	5
GB	Instructions for use	7
F	Mode d'emploi	10
I	Istruzioni per l'uso	13
P	Instrução de uso	16
E	Instrucciones para el uso	19

1

2

3

4

5

6

7


8

9

10

Cercon® PressMaster

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Directions for use

Transport and storage

Please store muffle and accessories away from the direct sunlight in a dry place.

Keep dry

Protect from sunlight

Please note the following symbols on the product labels:

REF Product number

LOT Batch number

Usable until

Obey instructions

Non-reusable

Indication

The Cercon® PressMaster universal pressing muffle allows, for the first time, fabrication of much longer span bridgework than is possible with other commercially available muffles for pressable ceramics. This enables bridges of up to 14 units to be pressed.

The unique design of the universal muffle allows the use of pellets of all geometries. The quadruple shamrock-like formation of the pressing sprues, similar to the cylinder of a revolver, allows the use of up to 8 pellets per 2.2 g in one single pressing. This is the equivalent of 6.4 grams of wax. In addition, you can use up to four different pellet shades simultaneously in one pressing. This puts pressable ceramics in a whole new perspective and offers an added ergonomic benefit for the technician through shortened work procedures.

Cercon® PressMaster is 90 mm in diameter and can be used in Cergo® press, Multimat® Touch and Press, Strato® Press (Elephant Dental), and EP 500 und 600 (Ivoclar) furnaces. Cercon® PressMaster is also suitable for use with zirconium oxide or galvanoformed frameworks.

Step by step instructions

considering Cercon® ceram press on Cercon framework as example

For detailed information please also consider the extensive directions for use of Cercon® ceram press.

General Notices

Always use a fireproof surface and never place a hot Cercon® PressMaster or other hot objects near flammable materials.

! IMPORTANT:

- Weigh the framework to determine the wax weight and ensure that fixing pins are used to prevent framework rotation.

1. Wax-Up (Fig. 1)

The wax-up must have a minimum layer thickness of 0.4 mm.

Before sprueing the restoration, weigh the pattern to determine the number of pellets needed.

The weight of the wax-up must not exceed 6.4g. This equals 17.6g of pressable ceramics (calculation: 0.8g wax = 2,2g pressable ceramic).

2. Sprueing (Fig. 2)

Attach the sprue to the restoration's vacuum form cap. This later guarantees easy removal of the muffle former.

Please note: Avoid sharp edges when attaching sprues! The pressing sprues should be no longer than 5 to 8 mm and at least 3.5mm in diameter.

2.1

Set the vacuum form cap with the sprued restoration on top of the muffle former (Fig. 3) and affix with soft wax. Firstly, apply a coat of vaseline to the walls of the muffle former and then place in the muffle form holder (Fig. 4).

2.2

Place the muffle former in the holder. Slide the sleeve over the muffle form holder (Fig. 5 and 6).

3. Investing

Investing requires 450 g of investment. Correct mixing of the investment is important; please observe the manufacturer's instructions.

After filling the muffle with investment (Fig. 7), place the spacer ring on top of the muffle while the investment is still fluid. This ensures a flat and stable muffle base for pressing (Fig. 8).

Remove the vacuum form cap after the investment has set.

Notice: This muffle has been tested explicitly with Cergo® fit Speed and Extra® Star in a practical rapid-heating environment.

4. Preheating

Please ensure the muffle is not placed with the base flush with the floor of the furnace as this can cause a steam pressure stasis during preheating. The large size of the Cercon® PressMaster requires a holding stage of 90 minutes at final temperature. The preheating temperature is defined by indication and material in use.

5. Pressing (Fig. 9)

To simplify placing several pellets, use the Fixpaste to glue them together. Each press sprue must be filled with the same number of pellets.

Notice: Directly before use, mix a small amount of powder with warm water on a separate glass plate or similar base.

The Cercon® PressMaster is compatible with the above mentioned furnaces. Centre the muffle on the left of the furnace! Please consider that Cercon® PressMaster holding and pressing times are twice the length of those of other press systems.

Pressing programme with

	Cergo® press:		Multimat® Touch & Press:	
	Cercon® ceram press	Cergo® Kiss	Cercon® ceram press	Cergo® Kiss
Starttemp.:	700 °C	700 °C	700 °C	700 °C
Vacuum:	cont	cont	50 hPa	50 hPa
Rate:	60 °C/min	60 °C/min	60 °C/min	60 °C/min
Final temp.:	940 °C	980 °C	940 °C	960 °C
Holding time:	40 min	40 min	40 min	40 min
Pressing:	20 min	40 min	20 min	40 min
Pressure:	5 bar	4.5 bar	3.0 bar	2.7 bar

The values indicated here are recommended values and serve only as a clue. Deviations of the burning results are possible. The burning results depend on the respective output per oven and are due to the manufacturer and to age. Therefore the recommended values must be individually adapted at each burning. We recommend a test burning to control the oven. All indications have been carefully elaborated and tested by us, but they are passed on without any guarantee.

6. Deflasking

Remove restorations from the muffle as normal. The Cercon® PressMaster has numbered segments for easier orientation.

Notes on the use of the universal muffle

Sprueing

Objects should be attached to the vacuum form cap with sprues 3.5 mm in diameter and not exceeding 5 to 8 mm in length. If longer sprues are required, as in longer span bridgework, the sprues should be interconnected to obtain a large flow area.

The inside of the muffle form holder has been provided with a coding that is repeated at regular intervals (1 to 4 dots) and replicated on the outside of the muffle to improve the pellet assignment process.

Pre-heating (Fig. 10)

The muffle should be placed on spacers in the pre-heating furnace, even if the floor of the heating chamber is corrugated (e.g. small open-web girders can be used for this), in order to allow the vapor pressure to deflect (especially if the speed investing technique was used).

Inserting pellets

Each sprue must be filled with the same number of pellets in order to prevent tilting of the investment plunger during the pressing process. This is essential even if objects were placed at only three sprues.

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DeguDent GmbH