DENTSPLY
TRUBYTE

Triad[®]

Visible Light Cure Provisional Material

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DIRECTIONS FOR USE

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CE
0086

CAUTION: U.S. Federal Law restricts this device to sale by, or on the order of, a licensed dental professional.

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DENTSPLY
TRUBYTE

Manufactured by:
DENTSPLY International Inc.
York, PA 17405-0872
1-800-786-0085

Authorized EU Representative
DENTSPLY DeTrey GmbH
Rodenbacher Chaussee 4
D-63457 Hanau-Wolfgang
Germany

www.trubyte.dentsply.com or www.dentsply.com

Indications for Use:

Triad Provisional Material is indicated for the fabrication of provisional restorations for teeth prepared for crowns, bridges, inlays, onlays and veneers.

Triad Provisional Material is a urethane dimethacrylate-based material that contains no methyl methacrylate monomer.

Contraindications:

Triad Provisional Material is contraindicated for patients and users with a history of allergic reaction to a) urethane resins; and/or b) methyl methacrylate monomer (contained in Triad VLC Bonding Agent only).

Warnings:

Triad Provisional Material and Triad VLC Bonding Agent contain polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice.

Precautions:

1. Users with special skin problems, cuts, or abrasions should wear protective gloves.
2. When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
3. All cure times cited are appropriate for Triad 2000 and Triad II curing units. Consult your Technique Manual or call Technical Support for Triad I cure times.
4. Store Triad Materials at or below 24°C/75°F and away from direct sunlight. Triad Materials may be refrigerated to extend shelf life.

5. Unused Triad Materials should be polymerized prior to disposal.
6. Do not use Triad Materials that crack or craze when handled.

Adverse Reactions:

1. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
2. Dust will be generated when grinding Triad Materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.
3. Some patients may experience an unpleasant taste from uncured Triad Materials.

STEP-BY-STEP INSTRUCTIONS

Laboratory Phase:

1. Create a preoperative cast from an accurate impression.
2. Fill edentulous spaces and temporarily restore fractured teeth so that the adjusted cast will represent the final restoration. Form a clear matrix over this preoperative cast. (See alternative C&B technique suggestions for options.)
3. Trim the transparent matrix to include the area being restored. It may be desirable to include adjacent areas to facilitate placement and maintain occlusion. Trim the matrix close to the gingival margins.

Clinical Phase:

1. Complete the tooth preparation and lubricate the teeth with petrolatum or saliva.
2. Select the desired shade of Triad Provisional Material and place it into the clear matrix in the areas to be restored. Place the filled matrix over the pre-

pared teeth. Press into place using the adjacent areas as a guide. Guide the patient into occlusion gently to enhance the seating process. Excess Triad Provisional Material will flow out around the margins.

In large, flat areas, cut a relief hole into the matrix to facilitate seating.

3. Lift the matrix partially off the tooth and return to the fully seated position three (3) or four (4) times. This step "burnishes" the surfaces of the resin material over the prepared tooth, helps eliminate undercuts, and allows better marginal adaptation.
4. Establish the fit, and remove excess material while the Provisional Material is soft. Then re-seat the matrix and use a hand-held visible light curing unit to initially set (harden) the provisional. (Interproximal excesses should be removed prior to light curing to facilitate removal and seating of the cured provisional.)
5. *Cure the facial and lingual surfaces with 30-second bursts from the hand-held light.* Remove the matrix and provisional. Evaluate the internal surface to assure that the material is hardened on the facial, occlusal and lingual surfaces, as well as on the surfaces in contact with the tooth preparation. The initial curing time with the hand-held light will depend upon thickness and length of the provisional being fabricated.
6. *Place the partially cured provisional into the Triad curing unit and cure for two (2) minutes.* Following this initial two (2) minute cure, remove the provisional restoration from the clear matrix and re-seat onto

the prepared teeth to check final occlusion and contact areas.

7. Trim with acrylic burs, and if necessary, relieve the internal surface slightly to allow space for the temporary cement.
8. When the occlusion, contact areas and margins are satisfactory, coat the provisional with Triad Air Barrier Coating (ABC) and *place in the Triad curing unit for a final cure of eight (8) minutes.*
9. Remove the provisional restoration and scrub gently with a soft brush and warm water to remove all traces of Triad ABC. Finish with polishing points, pumice or a rag wheel and polishing compound.
10. Place and cement the restoration as usual.

Instructions For Veneers:

Etch three small dots on the prepared enamel surface of the tooth to be restored (using usual procedures), and dry. Triad Provisional Material will adhere to these three dots sufficiently to provisionally restore the teeth until the next appointment, at which time the material may be removed. Note: Use of a bonding agent may make the Triad Provisional Material difficult to remove.

Place increments of the Triad Provisional Material onto the prepared surface, adapting it against the cut surface and sculpting the desired anatomy. Overbuild slightly to allow for finishing, as required. *Cure for 60 seconds using overlapping bursts from a hand-held light, just like for a direct-bonding restorative resin.*

Polish using normal office procedures for direct-fill resins.

NOTE:

1. Triad Provisional Material must be polished. When it cures, a thin, oxidized (air inhibited) layer forms on the surface. This needs to be polished away.
2. Do not use Triad ABC in the mouth.
3. Material may be added to the provisional restoration by using Triad VLC Bonding Agent. See the Triad Technique Manual for instructions on use of Triad VLC Bonding Agent.

INSTRUCTIONS FOR PROVISIONAL INLAYS AND ONLAYS

Inlays and onlays may be fabricated with Triad Provisional Material using a direct technique. Once the material is placed in the cavity preparation (following appropriate basing and lining of the preparation), interproximal areas should be checked, as in the crown provisional technique. This can help avoid polymerization of material in undercut areas.

Following placement and rough contouring of the material, use the hand-held curing light to harden the material (15 second overlapping bursts on each surface). Remove the material from the preparation, and then *cure in the Triad Curing Unit as per Steps 6 through 10, Clinical Phase.*

ALTERNATE PROVISIONAL CROWN/BRIDGE TECHNIQUES

As an alternative to fabricating a clear matrix on a cast, a thermally deformed plastic "button" or wax may be used directly in the mouth prior to tooth preparation. This alternative matrix can be lubricated with petrolatum and used as the matrix described in the clinical phase of the Triad Provisional technique.

Also, a visible light cure resin which is translucent (such as Triad Resiline) may be used to form the preopera-

tive matrix. Because these materials are composite materials, great care must be taken to lubricate the matrix with Triad Model Release Agent (MRA) so that the Triad material does not adhere to the matrix material.

DIRECT MOULDING

Triad Provisional Material also lends itself to direct moulding in the mouth. Remove sufficient Triad Provisional Material from the rope and place it on the prepared tooth. Use a clear mylar strip to assist in preventing material from being expressed into the interproximal areas. Direct carving can remove excess and begin to sculpt the desired tooth anatomy.

Establish the occlusion in this non-matrix technique by forming a functional general path directly in the uncured Triad Provisional Material. Be sure the opposing teeth are moist or lubricated. Moisture will not affect the cure of the Triad material.

When the contours and occlusion are satisfactory, harden the *Provisional Material with the hand-held light using 15-second multi-directional bursts.* Next remove the provisional, check for complete hardening, and remove flash as necessary. *Cure in the Triad Unit for 2 minutes.* Re-seat the provisional and make final adjustments with burs.

Finally, coat the provisional with Triad ABC and *cure it in the Triad Unit for eight (8) minutes.* Remove from the Triad Unit and wash away all traces of Triad ABC using warm water and a soft brush. Polish and cement.

INDIRECT PROVISIONAL

Methods similar to the clinical procedures may be used to fabricate provisionals *on a cast* rather than intra-orally. This gives the most predictable results for large/long-span restorations.

MATERIAL SAFETY DATA SHEET

FORM NO.: 85 SECTION I - GENERAL INFORMATION

1.1 Product Trade Name: Triad[®] VLC Provisional Material
 1.2 Part (Item) Number: 89302, 89303, 89304, 89305, 89306, 89307, 89308, 89309, 89310, 89311, 89349, 89732, 89733, 89734, 89735, 90130
 1.3 Chemical Name: N/A
 1.4 Synonyms: N/A
 1.5 Division Name: DENTSPLY Trubyte
 Address: 570 W. College Ave.
 City, State, Zip.: York, PA 17405-0872
 1.6 Phone No. for Info.: 717-845-7511
 Dentsply

"FOR CHEMICAL EMERGENCY"
 SPILL, LEAK, FIRE, EXPOSURE OR ACCIDENT
 CALL CHEMTREC - DAY OR NIGHT

1.7 Emergency Phone No.: 800-424-9300 1.12 National Stock No.: N/A
 Chemtrec 1.13 Specification No.: N/A
 1.8 Date: 11/20/96 1.14 Chemical Family: Acrylic
 1.9 Dentsply Contact: Scott E. Shaffer
 1.10 Signature: *Scott E. Shaffer* 1.15 Type/Grade/Class: N/A
 1.11 Formula: N/A 1.16 NRC/EPA Number: N/A

SECTION II - HAZARDOUS INGREDIENTS

2.1 CAS Number	Chemical Name (Ingredients)	%	TLV	OSHA PEL
NONE				

SECTION III - PHYSICAL PROPERTIES

3.1 Boiling Point: N/A.	3.11 Specific Gravity: 1.2
3.2 Vapor Pressure: <1 mm Hg	3.12 pH: N/A.
3.3 Vapor Density: >1	3.13 Evaporation Rate: N/A.
3.4 Critical Temp: N/A.	3.14 Viscosity: N/A.
3.5 Decomposition Temp: N/A.	3.15 % Volatile by Vol: Negligible
3.6 Melting/Freezing Pt.: <-40°F	3.16 Magnetism: N/A.
3.7 Solubility in Water: Negligible	3.17 Autoignition Temp: N/A.
3.8 Critical Pressure: N/A.	3.18 Corrosion Rate: N/A.
3.9 Perm. Exp. Limit: N/A.	3.19 Molecular Weight: N/A.
3.10 Appearance and Odor.: Tooth-colored, dough-like material in the uncured state. Hard plastic when cured.	

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

- 4.1 Flash Point (Method Used): 207°C (405°F) Cleveland Open Cup Method
 4.2 Flammable (Explosive) Limits:
 LEL: N/A. UEL: N/A.
 4.3 Extinguishing Media: Dry chemical, CO₂, Foam, Water Spray, Water Fog.
 4.4 Special Fire Fighting Procedures and Equipment: Fireman should wear protective clothing and self-contained respirators if large amounts are involved.
 4.5 Unusual Fire and Explosion Hazards: Irritating gases may be released.

FTriad® VLC Provisional Material, cont'd
SECTION V - HEALTH HAZARD DATA

- 5.1 Effects of Acute and Chronic Over Exposure: Material is nontoxic. Skin irritation is possible in sensitive individuals. Contact dermatitis (skin sensitization) can develop in sensitive individuals.
- 5.2 Emergency and First Aid Procedures: Eyes: Flush immediately with water for a minimum of 15 minutes and get medical attention. Skin Irritation: Wash with soap and water, remove contaminated clothing. Sensitization: Avoid further contact.

SECTION VI - REACTIVITY DATA

- 6.1 Unstable: Stable: X
- 6.2 Conditions to Avoid (Stability): High temperatures, direct sunlight.
- 6.3 Incompatibility (Materials to Avoid): Oxidizing agents.
- 6.4 Hazardous Decomposition Products: Acid fumes and smoke may result during initial stages of a fire.
- 6.5 Hazardous Polymerization: May Occur: Will Not Occur: X
- 6.6 Conditions to Avoid (Polymerization): High temperatures (polymerizes to form a hard plastic in 2 hours at 212°F (100°C) with a low exotherm).

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

- 7.1 Handling and Storage Precautions: Do not store near high temperature of ignition sources. Refrigeration prolongs shelf life. Store away from food and beverages.
- 7.2 Other Precautions: This material should be handled and used only in a manner consistent with the product instructions. Material is sensitive to visible light and will polymerize to form a hard plastic. Steps to be Taken in Case Material is Released or Spilled: No special procedures required; follow waste disposal procedure in 7.4.
- 7.3 Waste Disposal Procedures: Spread out materials to expose surfaces to sunlight or artificial light until solid, then dispose as a nonhazardous waste.
- 7.4 Waste Disposal Procedures: Spread out materials to expose surfaces to sunlight or artificial light until solid, then dispose as a nonhazardous waste.
- 7.5 Neutralizing Agent: N.A.

SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

- 8.1 Respiratory Protection: Use dust mask while grinding polymerized (cured) material to avoid breathing dust.
- 8.2 Ventilation: Use local exhaust when grinding polymerized (solid or cured) material.
- 8.3 Protective Gloves: Users with special skin problems, cuts, or abrasions should wear protective gloves (latex or PVC).
- 8.4 Eye Protection: Use safety glasses while grinding.
- 8.5 Other Personal Protective Equipment: None necessary.

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