**DENTSPLY**  
TRUBYTE

# Luci-Sof™

Soft Denture Liner System

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## DIRECTIONS FOR USE

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**CAUTION:** U.S. Federal Law restricts this device to sale by, or on the order of, a licensed dental professional.

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**DENTSPLY**  
TRUBYTE

Manufactured by:  
DENTSPLY Trubyte  
DENTSPLY International Inc.  
York, PA 17405-0872  
1-800-786-0085

Authorized EU Representative  
DENTSPLY DeTrey GmbH  
Rodenbacher Chaussee 4  
D-63457 Hanau-Wolfgang  
Germany

[www.trubyte.dentsply.com](http://www.trubyte.dentsply.com) or [www.dentsply.com](http://www.dentsply.com)

© 1997, 1999, 2002 DENTSPLY International Inc.  
All rights reserved.

Printed in U.S.A. 902860 Rev. F (8/02)

USA GB

**Indications for use:**

The Luci-Sof Soft Denture Liner System is intended for use as a long-term resilient lining material for removable dentures.

The Luci-Sof Denture Liner System consists of two components:

1. Luci-Sof Soft Denture Liner: a heat polymerizable, one component silicone putty that is based on dimethylpolysiloxane, dimethyl, methyl hydrogen siloxane copolymer; reinforced with amorphous silica. Benzoyl peroxide is the initiator for polymerization.
2. Luci-Sof Bonder: a heat polymerizable methacrylate monomer liquid and methacrylate polymer blend.

**Contraindications:**

1. Luci-Sof Bonder is contraindicated for patients and users with a history of allergic reactions to methyl methacrylate monomer.

**Warnings:**

1. Luci-Sof Bonder contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible individuals. Wash thoroughly with soap and water after contact. If skin sensitization or other allergic reactions occur, discontinue use. If dermatitis or other symptoms persist, seek medical advice.
2. Avoid inhalation or ingestion of

Luci-Sof Bonder. High vapor concentration can induce headache, nausea, and irritation of eyes and respiratory system. Liquid contact with eyes may cause possible corneal damage. Excessive long-term exposure may be associated with other more serious health effects. Monitor air quality per OSHA standards.

**Inhalation:** Move subject to fresh air. Give oxygen or artificial respiration as required.

**Ingestion:** Contact your regional Poison Control Center immediately.

**Eye Contact:** Flush eyes promptly with copious amounts of water for 15 minutes, and consult a physician. Wash skin with soap and water.

**Precautions:**

1. Do not use Luci-Sof on dentures that will result in denture base acrylic that is less than 2 mm in thickness.
2. Do not use the alternative long cure cycle (9 hours at 163° F) when lining a new denture.
3. Users with special skin problems, cuts or abrasions should wear protective gloves.
4. When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
5. Refrigeration of Luci-Sof Soft Denture Liner is recommended. Always bring to room temperature

before use. Do not store next to amine containing chemicals (ie. acrylic repair liquid, adhesive accelerator, etc.) as this may cause discoloration of the Luci-Sof material.

6. Store Luci-Sof Bonder below 80° F (27° C) and away from direct sunlight. May be refrigerated to extend shelf life.
7. Eugenol containing dental materials should not be used in conjunction with this product since they may interfere with hardening and bonding.
8. Unused Luci-Sof Soft Denture Liner should be polymerized prior to disposal.
9. Luci-Sof Bonder is a hazardous material. Dispose of in accordance with Federal, State, and Local regulations.


**Adverse Reactions:**

1. Corneal damage, headache, nausea, and vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precautions statements.)
2. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
3. Particulates will be generated when grinding prosthodontic materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

**STEP-BY-STEP INSTRUCTIONS:****Relining An Existing Denture**

Luci-Sof Soft Denture Liner System is a long-term silicone heat-cured soft denture liner material and bonder.

1. Pour a stone cast of the impression in the denture to achieve the desired contour.
2. Trim excess impression material, smooth with wax and finish to final shape around the periphery. **This preparation will help to minimize finishing later.**
3. Invest.
4. Open flask and remove wax with boiling water.
5. Remove impression material.
6. Wash with detergent and water and rinse with clean boiling water.
7. Grind the denture surface with a clean carbide cutter to remove tray adhesive and sufficient acrylic. **A clean surface is necessary to achieve a superior bond to the denture base.** Avoid knife edges. Grind flange border to a 90° butt joint.

8. Shake the Luci-Sof Reline Bonder well and generously coat the areas

Luci-Sof™ Soft Denture Liner, cont'd  
SECTION V - HEALTH HAZARD DATA

- 5.1 Effects of Acute and Chronic Over Exposure: None known.
- 5.2 Emergency and First Aid Procedures: INHALATION: No irritation to eyes or respiratory passages. No injury is likely from relatively short exposure of less than 8 hours. No first aid should be needed. EYES: Direct contact may cause temporary discomfort with mild redness similar to windburn. Immediately flush with water. SKIN: A single prolonged exposure (24 to 48 hours) caused no known adverse affects. INGESTION: Small amounts transferred to the mouth by fingers during use, etc. should not injure. Swallowing large amounts may cause digestive discomfort. No first aid should be needed.

SECTION VI - REACTIVITY DATA

- 6.1 Unstable: Stable: X
- 6.2 Conditions to Avoid (Stability): Water, Alcohols, acidic or basic materials, and many metals or metallic compounds, when in contact with product, liberate flammable hydrogen gas, which can form explosive mixtures in the air.
- 6.3 Incompatibility (Materials to Avoid): Oxidizing materials can cause a reaction. See conditions to avoid.
- 6.4 Hazardous Decomposition Products: Silicon Dioxide, Carbon Dioxide, and traces of incompletely burned carbon products.
- 6.5 Hazardous Polymerization: May Occur: Will Not Occur: X
- 6.6 Conditions to Avoid (Polymerization): N.A.

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

- 7.1 Handling and Storage Precautions: Store in cool dry place.
- 7.2 Other Precautions: Store away from food and beverages.
- 7.3 Steps to be Taken in Case Material is Released or Spilled: Collect and contain for salvage or disposal.
- 7.4 Waste Disposal Procedures: Follow Federal, State and Local regulations for disposal.
- 7.5 Neutralizing Agent: N.A.

SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

- 8.1 Respiratory Protection: None should be needed.
- 8.2 Ventilation: Local exhaust not needed.
- 8.3 Protective Gloves: For users with skin problems, cuts or abrasions - wear protective gloves.
- 8.4 Eye Protection: Use safety glasses while grinding.
- 8.5 Other Personal Protective Equipment: None

SECTION IX - TRANSPORTATION

- 9.1 Applicable Regulations... IMOC: N.A. IATA: N.A. 49 CFR: N.A.  
Military Air (AFR 71-4): N.A. Tariff 6D: N.A.
- 9.2 Shipping Name: N.A.
- 9.3 I.D. Number: N.A. 9.9 Unit Container: N.A.
- 9.4 Report Quantity: N.A. 9.10 DOT SPSC Container: N.A.
- 9.5 Hazard Class: N.A. 9.11 DOT Exempt/DOD/CCN: N.A.
- 9.6 Labels: N.A. 9.12 Limited Quantity: N.A.
- 9.7 Aerosol Propellant(s): None
- 9.8 Net Explosive Weight: N.A.

MATERIAL SAFETY DATA SHEET

FORM NO.: 180 SECTION I - GENERAL INFORMATION

- 1.1 Product Trade Name: Luci-Sof™ Bonder
- 1.2 Part (Item) Number: 904556, 904559
- 1.3 Chemical Name: Methyl Methacrylate Monomer
- 1.4 Synonyms: MMA
- 1.5 Division Name: DENTISPLY Trubyte
- Address: 570 W. College Ave.
- City, State, Zip: York, PA 17405-0872
- 1.6 Phone No. for Info: 717-845-7511

"FOR CHEMICAL EMERGENCY"

SPILL, LEAK, FIRE, EXPOSURE OR ACCIDENT

CALL CHEMTREC - DAY OR NIGHT

- 1.7 Emergency Phone No.: 800-424-9300 1.12 National Stock No.: N.A.
- Chemtrec 1.13 Specification No.: N.A.
- 1.8 Date: 8/13/02 1.14 Chemical Family: N.A.
- 1.9 Dentsply Contact: Scott E. Shaffer Methacrylic ester
- 1.10 Signature: *Scott E. Shaffer* 1.15 Type/Grade/Class: N.A.
- 1.16 NRC/EPA Number: N.A.
- 1.11 Formula: CH<sub>2</sub>(CH<sub>3</sub>)COOCH<sub>3</sub>

SECTION II - HAZARDOUS INGREDIENTS

CAS Number	Chemical Name (Ingredients)	%	TLV	OSHA PEL
80-62-5	Methyl Methacrylate	<79	100ppm	100ppm
97-90-5	Ethylene Glycol Dimethacrylate	<6.0	1mg/m <sup>3</sup>	
94-36-0	Benzoyl Peroxide	<0.5	5mg/m <sup>3</sup>	

SECTION III - PHYSICAL PROPERTIES

- 3.1 Boiling Point: 214°F/101°C 3.11 Specific Gravity: 0.95 (820°C/60°F)
- 3.2 Vapor Pressure: 29mmHg@68°F 3.12 pH: N.E.
- 3.3 Vapor Density: 3.45 Air=1 3.13 Evaporation Rate(BAC=1): >1 Bac=1
- 3.4 Critical Temp.: N.A. 3.14 Viscosity: N.A.
- 3.5 Decomposition Temp.: N.E. 3.15 % Volatile by Vol.: 85%
- 3.6 Melting/Freezing Pt.: N.A. 3.16 Magnetism: N.A.
- 3.7 Solubility in Water: Slight 3.17 Autoignition Temp.: 815°F(435°C)
- 3.8 Critical Pressure: N.A. 3.18 Corrosion Rate: N.A.
- 3.9 Para. Exp. Limit: N.A. 3.19 Molecular Weight: N.A.
- 3.10 Appearance and Odor: Cloudy viscous liquid, characteristic acrylic odor.

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

- 4.1 Flash Point(Method Used): 55°F(13°C) TOC, 55°F(13°C) TCC
- 4.2 Flammable (Explosive) Limits:  
LEL: 2.12% UEL: 12.5%
- 4.3 Extinguishing Media: CO<sub>2</sub>, foam, water spray, and water fog.
- 4.4 Special Fire Fighting Procedures and Equipment: Water may be ineffective unless used as a fine spray or fog. Use water spray to cool fire-exposed containers of methyl methacrylate. Self contained breathing apparatus and protective clothing should be worn to fight fire if large amounts of liquid are combusting.
- 4.5 Unusual Fire and Explosion Hazards: Vapors may travel to source of ignition and flash back. Avoid ignition sources or excessive temperatures. Heat can induce polymerization with rapid release of energy and closed containers may rupture explosively (avoid by cooling containers with water spray).

Luci-Sof™ Bonder, cont'd  
SECTION V - HEALTH HAZARD DATA

- 5.1 **Effects of Acute and Chronic Over Exposure:** Inhalation: High vapor concentrations can induce headache, nausea, and irritation of respiratory system. Eye Contact: Liquid contact with eyes will cause irritation and possible corneal damage. Skin Contact: Material may cause skin irritation. Can cause contact dermatitis in some individuals. Avoid contact if allergy/dermatitis develops.
- 5.2 **Emergency and First Aid Procedures:** INHALATION: Move subject to fresh air. INGESTION: If conscious, give 2 glasses of water and consult a physician immediately. Give oxygen or artificial respiration as required. EYE CONTACT: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician if irritation is apparent. SKIN CONTACT: Wash skin with soap and water. Contact physician if irritation persists.

SECTION VI - REACTIVITY DATA

- 6.1 **Unstable:** Stable: X (when stored properly)
- 6.2 **Conditions to Avoid (Stability):** Elevated temperatures and direct sunlight.
- 6.3 **Incompatibility (Materials to Avoid):** Reducing or oxidizing materials, acids, bases.
- 6.4 **Hazardous Decomposition Products:** None known.
- 6.5 **Hazardous Polymerization:** May Occur: X Will Not Occur:
- 6.6 **Conditions to Avoid (Polymerization):** Elevated temperatures, oxygen-free atmosphere inhibitor depletion (due to excessive aging), direct sunlight, contamination with polymerization catalysts.

SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

- 7.1 **Handling and Storage Precautions:** Store in cool dry place out of direct sunlight. Use adequate ventilation.
- 7.2 **Other Precautions:** Avoid repeated breathing and contact with the skin. Keep containers closed. Store away from food and beverages. Keep away from sparks and open flame.
- 7.3 **Steps to be Taken in Case Material is Released or Spilled:** Eliminate ignition sources. Avoid breathing vapors. Dike & absorb spill with inert material (sand, soda, ash, vermiculite, etc.) and transfer to containers for disposal; use non-sparking tools.
- 7.4 **Waste Disposal Procedures:** After addition of excess inhibitor, incinerate under controlled conditions in a safe open area, according to federal, state and local regulations
- 7.5 **Neutralizing Agent:** N.A.

SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

- 8.1 **Respiratory Protection:** Self contained breathing apparatus during high vapor concentration.
- 8.2 **Ventilation:** Local exhaust to keep vapor below TLV and LEL. Mechanical ventilation should have explosion proof motors.
- 8.3 **Protective Gloves:** Impervious gloves (butyl rubber preferred).
- 8.4 **Eye Protection:** Splash proof goggles.
- 8.5 **Other Personal Protective Equipment:** Wear impervious apron and overshoes during clean-up operation.

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