## DESCRIPTION

The Radica Provisional System includes visible light-cured, resin-based composites and accessories used to create provisional restorations for patients awaiting final prosthodontic restoration.

## INDICATIONS

Indicated for fabrication of temporary restorations

## CONTRAINDICATIONS

1. Radica is contraindicated for patients with a known sensitivity to acrylic resins.
2. Radica is contraindicated for fabricating permanent restorations.

## WARNINGS

1. Radica contains components that can cause irritation to the eyes and skin. Prolonged contact may cause skin sensitization (allergic contact dermatitis). On direct contact with the skin, wash thoroughly with soap and water. Consult a doctor if dermatitis or other symptoms persist.
2. The sealer is highly flammable and can cause irritation to the eyes and skin. The sealer can cause serious injury if it comes in contact with the eyes. Irritant to respiratory system and skin. May cause allergic skin reactions. May cause sensitization on repeated skin contact particularly with susceptible people. Use only in a well-ventilated workspace to avoid inhalation. High concentrations of vapours from sealer/stains in a poorly ventilated area may induce headaches, dizziness and irritation to the eyes and respiratory system. Prolonged excessive exposure can cause other serious health problems. Monitor the air quality according to OSHA standards.

First aid in case of:

- *Eye contact:* Rinse the eyes immediately with plenty of water for at least 15 minutes; consult a doctor.
- *Inhalation:* Move the patient into the fresh air. Give oxygen or artificial respiration if necessary.
- *Ingestion:* If fully conscious, the patient should rinse his/her mouth and drink 2–4 cups of milk or water. Do not induce vomiting. Seek medical attention.
- *Skin contact:* Wash the skin for at least 15 minutes with soap and water. Seek medical attention if symptoms persist.

## PRECAUTIONS

1. Only use the Radica system in a well-ventilated workspace.
2. This product is light sensitive! Protect from light and always close the container to ensure that the material properties are not impaired.
3. Store the product in a cool, dry place at 15–26°C.
4. Keep sealer bottles tightly closed when not in use. The material is volatile and flammable. Keep away from heat, sparks and naked flames.
5. Avoid inhalation of grinding dust; use an extractor and an industrial ventilation system.
6. Restorations with a mesial to distal span greater than 38 mm cannot be polymerized in the Triad® 2000 light-curing unit. Restoration/model with a height greater than 40 mm should not be processed in the Eclipse Junior VLC Curing Unit.
7. When repairs or adjustments are to be made intraorally by the dentist, ensure that the cured sealer layer is completely removed. Repolish using a manual polishing technique after repair.

## SIDE EFFECTS

The Radica system contains components that may cause dry skin, sensitization (allergic contact dermatitis) or other allergic reaction in susceptible persons.


## EQUIPMENT

### Radica syringe warmer

Radica resin must be heated prior to application. Connect the syringe warmer and switch it on. The warmer is preset to 60–64°C. Allow the warmer to heat up for approx. 30 minutes. Place the required dentine and enamel syringes in the holders of the warmer approx. 10 minutes prior to use. The Radica syringes can remain in the warmer during the day, but should not be kept heated for an extended period of time.

### Electric wax knife

The Ney® electric wax knife is used for contouring the resin. The special nickel-coated tips prevent discoloration of the resin. Always use electric wax knives with chromium or nickel-coated steel tips. Keep the wax tips clean. Do not use hot instruments that have been heated by an external source.

### Light-curing units


Only the Eclipse Junior VLC Curing Unit, the Triad® 2000 light-curing unit (maximum span of 38 mm with the Triad 2000) or the Eclipse® polymerization unit should be used.

Adhere to the following guidelines to ensure optimum curing.

1. When using the Eclipse Junior VLC Curing Unit, the combined height of the model and restoration should not exceed 40 mm. The maximum height is indicated by a perforated line on the Eclipse Junior VLC Curing Unit.
2. When using the Triad 2000 VLC unit, the RESTORATION, not the MODEL, must be placed in the centre of the platform below the lamp housing at the height indicated by the red line. See the photographs below.

*Eclipse Junior VLC Curing Unit*

*Triad 2000*

*Triad 2000  
Note: centred and at the  
correct height for curing*

## PREPARING THE MATRIX

1. Ensure that the teeth to be prepared are clean and free of debris.
2. Close any edentulous areas with denture teeth or a wax-up.
3. Use Radica matrix putty; thoroughly mix the base and catalyst in a ratio of 1:1.
4. Apply the matrix putty to the area of the model to be used for the temporary restoration. Include the adjacent teeth and apply the material at least 2 mm over the marginal area of the pattern.
5. Allow the material to cure. Remove the matrix from the model.
6. Cut the matrix with a sharp instrument to within approx. 2 mm of the area to be restored.

*Closing the edentulous space*

*Matrix on the model*


*Cutting the matrix*

## DIE /MODEL PREPARATION

*Note: It is important to reduce the dies correctly. If the dies are reduced too much, the temporary will be too thick and will have to be adjusted by the operator prior to fitting. The minimum thickness of connectors should be 4x5 mm.*

1. Use a Midwest® HP 1558 carbide bur with a 1.0 mm diameter. Prepare depth cuts in the model stone as a guide.
2. Use a model saw or diamond disc to open up the contact areas to the adjacent teeth by 1 mm. Ensure that the adjacent teeth are not damaged. Reduce the surfaces by approx. 1.0 mm with a carbide bur.

*Midwest 1558 bur with a 1.0 mm diameter*

*Preparing depth cuts*

*Opening up the interproximal area*

*Surfaces (crown walls) reduced by 1.0 mm*

*Applying separating agent*

## **RADICA™-DENTINE & ENAMEL SHADE COMBINATION CHART**

Dentine	A1	A2	A3	A3.5	B1	C2	C4	D2	i2
Enamel	Light	Light	Light	Light	Light	Light	Medium	Light	Light


### **FABRICATING THE TEMPORARY**

1. To build up the incisal, syringe a small amount of preheated enamel resin into the matrix. Contour the resin with the electric wax knife and remove any excess material. Allow the resin to air cool for approx. 1 minute.
2. The 3 dentine shade modifiers (red-orange, yellow-orange and pink) are used for customizing the restoration and can also be applied in the matrix.
3. Then quickly syringe dentine into the matrix. Place the matrix immediately onto the prepared model; any excess material exudes from the margins.
4. Allow the material to cool for approx. 2–3 minutes.
5. Remove the matrix from the model; if necessary enamel or dentine can be added to the restoration using the electric wax knife.
6. Excess resin is easily removed with an instrument. The connectors should be U-shaped (not V-shaped or tapered) and have a minimum diameter of 4 x 5 mm.
7. Apply a thin film of sealer with a brush. Ensure that the sealer does not overlap onto the model when applying it to the restoration, as this would make the restoration difficult to remove from the model following polymerization.

*Note: Clean the brush with acetone immediately after use; if not ensure it is protected from the light.*

*Applying the enamel*

*Contouring the enamel*

*Applying the dentine modifier*

*Completing with dentine*

*Repositioning the matrix*

*Removing the matrix*

8. Place the model and restoration in the Eclipse Junior, Triad 2000 or Eclipse polymerization unit according to the instructions for the unit and complete an initial cure cycle according to Table 1.

*Adding to the pattern*


*Removing excess*

*Excess resin removed*

*Applying the sealer*

*Curing the resin*

9. When fabricating crowns, proceed to FINISHING THE TEMPORARY RESTORATION, STEP 1.
10. With bridges use an instrument to remove the temporary restoration from the model following polymerization. For dark shades, e.g. A3.5 or C4, the bridge should be cooled before it is removed from the model (e.g. with compressed air).
11. Apply a thin film of sealer to the underside of the bridge pontic with a brush.
12. Place the bridge upside down on the turntable of the Eclipse Junior or Eclipse polymerization unit; with the Triad 2000 place it in the centre of the turntable and start the Pontic cure programme. Adhere to the times at CURING BRIDGES in Table 1.

*Removing the bridge  
from the model*

*Applying the sealer*

## **FINISHING THE TEMPORARY RESTORATION**

1. Remove the restoration from the model with an instrument. Any additional adjustments to the surface can now be made.
2. Clean the restoration by sandblasting all the surfaces (inner and outer) with aluminium oxide at 1.38 bar (20 psi). Then clean all the surfaces with a steam cleaner and allow to dry.
3. Apply a thin film of sealer to the outer surfaces of the restoration with a brush. (Clean the brush immediately with acetone.)  
If the bridge is to be stained, apply the stains immediately using a fine brush.
4. There are five stains (brown, blue, white, orange and olive-green); these can be used for surface characterization.

**Technique:** Place a small amount of stain on a mixing block, mix it with a small amount of sealer and then apply the stain.


5. Cure the sealer in the Eclipse Junior, Eclipse or Triad 2000 polymerization unit.  
Adhere to the times at CURING THE SEALER in Table 1.
6. Radica restorations can also be polished mechanically.

*Sandblasting with aluminium oxide*

*Mixing the stains*

*Applying the stains*

*Applying the sealer*

*The finished temporary bridge*

## REPAIRS/ADD-ONS/TISSUE TINT

1. Roughen the area intended for repair, add-on or tissue tint using a carbide bur.  
Extend the roughened area approx. 1 mm beyond this area.
2. If necessary, apply separating agent to the relevant area of the model and then replace the restoration on the model.
3. Moisten the roughened surface with a thin film of sealer.
4. Dispense a small amount of dentine or enamel as required on a mixing block.  
Then apply the resin with an electric wax knife to the roughened surface.
5. Allow the resin to cool for 2 minutes.
6. Apply a thin film of sealer to the adjusted area using a brush.  
Clean the brush immediately with acetone; if not ensure it is protected from the light.
7. Cure in the Eclipse Junior, Triad 2000 or Eclipse polymerization unit.  
Adhere to the times at INITIAL CURE in Table 1.
8. Finish using the procedure described above.

*Roughening the surface*

*Applying the sealer*

*Adding to the restoration*

## RADICA DIAGNOSTIC SYSTEM

Radica can be used both for temporary restorations and diagnostic restorations. It produces a very aesthetic diagnostic restoration.

1. Reduce the teeth uniformly by 1.0 mm.
2. Wax up the pattern fully and optimize the dentition.
3. Take an impression of the pattern with putty.
4. Fabricate a restoration using Radica as described above.  
Follow the steps in the section "Fabricating the temporary".
5. Cure the diagnostic restoration according to instructions.

*Working model*

*Preparing depth cuts*

*Uniform 1.0 mm reduction*

*Wax pattern*

*Matrix*

*Finished diagnostic temporary restoration*

## BEFORE DELIVERY OF A RADICA RESTORATION TO THE PRACTICE

Ensure that the inner surfaces of the crowns are clean: sandblast them if required with aluminium oxide at 1.38 bar (20 psi).

*Before cleaning*

*After cleaning*

## TABLE 1. POLYMERISATION TIMES FOR DENTSPLY LIGHT-CURING UNITS

	Eclipse Junior VLC Curing Unit	Triad® 2000 light-curing unit	Eclipse® polymerization unit
Initial cure	5 min.	10 min.	3 min.
Curing bridges	1.5 min.	2 min.	1/2 min.
Curing the VLC Sealer	2 min.	2 min.	1 min.

*Note: Do not use early versions of the Triad light-curing unit.*

# OPERATING INSTRUCTIONS FOR THE ECLIPSE JUNIOR VLC CURING UNIT

## INSTALLATION

The display window lights up when the unit is connected.  
The version number of the software is displayed.  
The unit performs a self-test before it is ready for operation.  
The name "Eclipse" is displayed automatically.

## CHOOSING THE PROGRAMME

Press the F3 button at the right below the display – this button allows you to scroll through the available programmes.

**When "Radica" is displayed,** press the middle button F2. This switches the menu to the Radica programme.

## PROGRAMME SELECTION

The two outer buttons can be used for scrolling through the Radica programme selection.

Select **the "5:00" minute programme** for initial curing.

Place the restoration in the unit and close the door manually.

Press the "Start" button.

Select **the "1.5" minute programme** for curing pontics.

Select **the "2.00" minute programme** for curing the sealer and stains.

## TIPS, TRICKS & TROUBLESHOOTING

<i>Can a hand-held polymerization unit be used with Radica?</i>	<ul style="list-style-type: none"> <li>• Yes – but only to “tack” the resin in the correct position. A full curing cycle still has to be carried out.</li> </ul>
<i>The cured temporary restoration cannot be removed from the model.</i>	<ul style="list-style-type: none"> <li>• Insufficient separating agent has been applied to the model.</li> <li>• Undercuts on the model</li> <li>• Sealer has been inadvertently applied to the model.</li> </ul>
<i>Grey streaks/discoloration in the finished restoration</i>	<ul style="list-style-type: none"> <li>• An electric wax knife was used that did not have nickel or chromium-plated tips.</li> <li>• A warm instrument was used that was heated by an external source.</li> </ul>
<i>Cavities or porosity in the cured material</i>	<ul style="list-style-type: none"> <li>• Air was trapped in the resin when it was syringed into the matrix.</li> <li>• The surface was not prepared correctly when adding on.</li> </ul>
<i>Inaccurate shade reproduction</i>	<ul style="list-style-type: none"> <li>• Too much or too little enamel was used in proportion to the dentine.</li> <li>• Incorrect enamel was used.</li> <li>• Too much dentine modifier was used.</li> <li>• Incorrect positioning of the model during polymerization in the light-curing unit (only Triad 2000)</li> <li>• An approved light-curing unit was not used.</li> <li>• The dentist did not use the correct resin liner.</li> </ul>
<i>The resin cannot be dispensed from the syringe.</i>	<ul style="list-style-type: none"> <li>• The syringe was not left long enough in the warmer.</li> <li>• The syringe was not fully inserted into the warmer.</li> <li>• The tip of the syringe was not closed properly and the resin hardened.</li> </ul>
<i>The finished temporary restoration does not have the required strength.</i>	<ul style="list-style-type: none"> <li>• The correct curing cycle was not used.</li> <li>• A recommended light-curing unit was not used.</li> <li>• The recommended connector thickness was not used.</li> <li>• Incorrect positioning of the model in the light-curing unit.</li> <li>• The walls of the restoration are too thin.</li> </ul>
<i>The surface of the temporary restoration is rough after curing.</i>	<ul style="list-style-type: none"> <li>• The VLC sealer was already slightly cured before it was applied to the surface.</li> <li>• The surface was not smooth before the sealer/ glaze was applied.</li> <li>• The application brush was contaminated.</li> <li>• Brush strokes were applied to the sealer after it began to dry.</li> </ul>
<i>There is no contouring on the surface.</i>	<ul style="list-style-type: none"> <li>• The sealer was applied too thickly.</li> </ul>
<i>Blue discoloration on the surface</i>	<ul style="list-style-type: none"> <li>• Incorrect mixing ratio of putty base and catalyst</li> </ul>

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