# VERBLENDKUNSTSTOFFE

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## Instructions for use

### **1. Indications for use:**

Light-curing plastic for making individual impression spoons.

### MATERIAL

Triad® VLC tray/base plates are available in pink and neutral and consist of a matrix of multifunctional acrylates and inorganic fillers.

### PROPERTIES

Triad® VLC tray/base plates exhibit very

- minimal shrinkage,
- minimal water uptake and
- high dimensional stability.

### **2. Contraindication:**

Do not use this product for preventive reasons if patient has shown hypersensitivity to any of the components it contains.

### **3. Warnings:**

For dental use only.

### **4. Precautions:**

Please note:

Use a dust extractor.

Avoid contact with skin.

Wear suitable gloves.

### STORAGE

Do not store at temperatures exceeding 25 °C.

Higher temperatures affect the handling properties.

## 5. Adverse reactions:

Do not use the product in case of manifestations of irritation or known allergy to a component of the product.

## 6. Step-by-Step instructions:

1. Block out undercuts with wax. To fabricate a tray, place a sheet of wax on the model to create a spacer. Extend the wax to within 2 – 3 mm of the tray periphery. Stops may be cut into the wax on the alveolar ridges.
2. Apply either Triad® Model Release Agent (MRA) or an alginate separating agent such as Isolant.
3. Unpack the Triad® VLC tray/base plate, adapt it to the model carefully and trim the peripheries. The excess material is easily used to form a tray handle.
4. Coat the entire tray or base plate with Triad® ABC varnish.
5. Polymerize the tray or base plate on the model in a light curing unit. Polymerization time in a Triad® light curing unit: 5 minutes.
6. Release the tray or base plate from the model, coat the underside with Triad® ABC varnish and polymerize again for 5 minutes.
7. Brush off the Triad® ABC varnish under running water, dry the tray or base plate and trim as usual with tungsten carbide cutters and sandpaper.

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DeguDent GmbH