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Contex[®]

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

XI

Reizend
Irritant
Irritante
Irritante

DeguDent

A Dentsply International Company

Context®

D

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Context®

Context is an anti-flux for precious metal soldering. Precious metal areas covered with Context will not be affected by the solder. This avoids unintentional spatter of the solder onto areas close to the soldering seam.

Because of this, Context proves particularly suitable for the soldering of spring- and friction pins, anchor- and root pins, precision attachments and in all cases where precision soldering is necessary.

Context is for use in goldsmith work in conjunction with all red gold alloys and with palladium-bearing white gold alloys, but not with nickel-bearing white gold alloys.

Availability: Jar of 50 ml. REF 53 2524 0004

Instructions for use:

1. Cover precious metal areas which need protection against spreading of solder with Context by applying a medium layer with a thin instrument.
2. Heat object slightly (for instance with a small flame) until a dry, adhesive layer has formed.
3. If necessary, place object to be soldered in soldering investment.
4. Then apply flux (f.e. Oxyxon® or Anoxan®) onto the soldering gap which was kept free of Context.
5. Solder in the usual way.
6. Close Context container tightly after use and prevent from drying and frost.

Attention!

Irritating to eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water.

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DeguDent GmbH