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Deguvest[®] California

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

DeguDent

A Dentsply International Company

Deguvest[®] California

D	Gebrauchsanweisung	3
GB	Instructions for use	7
F	Mode d'emploi	11
I	Istruzioni per l'uso	15
E	Instrucciones para el uso	19

Deguvest[®] California

D

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Deguvest[®] California

– plaster-bonded investment compound for precision casting –

Deguvest California is a high-quality plaster-bonded, graphite-free investment compound. The fine grain size guarantees an outstanding surface quality with excellent model detail reproduction.

The product is also notable for the ease with which it can be removed from the casting ring, one of the most important considerations when judging an investment compound e.g. for inlay technique.

Deguvest California is suitable exclusively for low-melting precious metal dental alloys, e.g. Degunorm[®], Degulor[®] M and thus offers a useful addition to the GoldenGate[®] system.

When heating, the pre-heating temperature for the alloy being used must be taken into consideration whereas 750 °C should not be exceeded. Deguvest California is mixed with distilled water.

1 Technical specifications

Mixing ratio powder : distilled water	100 g : 32 – 40 ml
Total expansion (in % linear)	1.3 – 1.6
Processing range (min. at 20 °C)	7–8

2 Important notes

To ensure the high quality of Deguvest California, the bag should be well sealed after the powder has been removed.

It is recommended to fill the contents of the bag into an air tight container to protect the investment compound from humidity.

Avoid contact with phosphate-containing investment compounds as they could be damage as a result. For mixing purposes a separate bowl and spatula should therefore be used.

Store in a cool and dry place.

Deguvest California can be stored for 18 months.

- Fix on the sprue former with wax.
Ensure that there is approx. 7 mm distance between the model and the edge of the ring in order to stabilize the investment compound “bottom”.
Also ensure that no sharp edges are at the sprue junctions, otherwise the investment compound could be torn off by in-flowing molten alloy.
- Brush or spray wax modelation with Waxit® wetting agent and remove excess wetting agent (do not use compressed air!).
- Shake Deguvest California briefly in bag or stock container, remove and weigh precisely to the gram.
Pour the distilled water in accordance with the above mixing ratio into a clean and dry mixing bowl and sprinkle the powder.

Mix the investment compound by hand for approx. 30 seconds with a spatula, allow to stand under vacuum for approx. 30 seconds, then mix under vacuum for 60 seconds with a mixing and evacuation unit (Multivac®). 130 g Deguvest California are required to fill a 3x casting flask.

4.1 Investment

After mixing, place the mixing bowl on a vibrator (approx. 30 seconds) and allow the investment compound to flow together. Then let the compound to flow into the mould, shaking gently and ensuring that no bubbles form. After the mould has been filled, switch off the vibrator immediately.

4.2 Setting

Allow the mould to set for 45–60 minutes, depending on the size (see “Heat treatment” table, point 5.1). If the mould is preheated too soon the investment compound can crack.

5 Pre-heating /casting

5.1 Pre-heating

- After the recommended drying time the wax can be eliminated from the mould in accordance to the table below and heated to the pre-heating temperature recommended for the alloy being used. The maximum heat-up rate should not exceed 7 °C per minute.

Heating	Size of casting flask			
	1x	3x	6x	9x
Setting/drying time	45 min	50 min	55 min	60 min
Initial temperature of pre-heating furnace	200 °C	200 °C	200 °C	200 °C
Wax-elimination time (at 290 °C)	30 min	40 min	50 min	60 min
Holding time at final temperature	20 min	30 min	45 min	60 min
Ensure that the <u>maximum</u> pre-heating temperature is not higher than 750 °C!				

- Place the mould so into the furnace that the wax can easily run out. Wax residues could penetrate into the capillaries of the investment compound and carbonize.
The moisture of the investment compound can also be used to eliminate the wax. Moulds which have dried overnight should be soaked slightly in water before the casting rings are placed in the furnace to eliminate the wax.

5.2 Casting

Follow the instructions of the alloy manufacturer.

5.3 Devestment

Allow the casting ring to cool slowly off at room temperature and remove the investment compound under running water with a small brush (old toothbrush).

6 Safety information

Mixing with other liquids or materials health hazards as well as impairments to the cast quality cannot be excluded.

All dental investment compounds contain parts of crystalline quartz and/or cristobalite.

These substances can adversely affect health if they are breathed in. Always use a sufficiently powerful extractor or wear a mask.

Carry out the “coarse” removal of the casting under water or soak the investment material mould.

6.1 Information for the processing of mixing liquids for dental investment compounds

Following skin contact:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

For further details please refer to the safety data sheets.

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DeguDent GmbH