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Deguvest[®] SR

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso
Инструкция по
применению

DeguDent

A Dentsply International Company

Deguvest® SR

D	Gebrauchsanweisung	4
GB	Instructions for use	16
F	Mode d'emploi	29
I	Istruzioni per l'uso	42
E	Instrucciones para el uso	55
RUS	Инструкция по применению	68

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1	Before starting work	18
2	Technical information	18
3	Investment storage	19
4	Investment instructions	19
4.1	Preparation	19
4.2	A few words about measuring out liquid	20
4.3	Mixing powder and liquid	21
4.4	Investment	21
5	Stetting times	22
6	Preheating and casting	22
6.1	Conventional Heating	22
6.2	Quick Cast method	22
7	Casting	23
8	Devesting	23
9	Recommendations for Non-precious alloys	24
	Investing with Deguvest SR - Quick reference guide	25
	Tables	26

Deguvest® SR

"better results, healthy environment!"

CONGRATULATIONS!

In choosing the DUST REDUCED **Deguvest SR** you have opted for a carbon-free, phosphate-bonded and gypsum-free universal investment, which is suited to all types of precious metal casting.

WHY CHOOSE A DUST REDUCED INVESTMENT?

Respirable crystalline silica (which is found in every investment) is considered by most health and safety authorities to be carcinogenic and is but one substance on a long list of carcinogenic substances. We realised the dangers involved in working with such materials and therefore decided to develop a dust reduced investment.

Dental technicians are exposed to hazardous substances on a daily basis. Being under pressure in the laboratory causes us to neglect health and safety precautions. As we care about your health we at DeguDent engaged very intensively in the subject of dust prevention. **Deguvest SR** tackles this problem at it's source! Very low dust levels during handling means one less hazard for the technician at work!

QUALITY FROM START TO FINISH!

We employ a Quality Management System which enables us to control and define the quality of each product at each station during the production process.

The finished product is placed under extreme scrutiny before finally being passed for use.

This gives the technician in the field the security that our product will give the usual top results – time after time – and – charge after charge!

APPLICATION

Deguvest SR can be heated in the *conventional manner*, using holding stages and then heating to final temperature (casting temperature) or, alternatively, using the "*Speed-Cast*" method i.e., placing the casting ring, directly in the furnace after setting, at final casting temperature.

1 Before starting work ...

we would like you to read through the following safety precautions.

Even though **Deguvest SR** has an extremely low dust development, we would still recommend (in the interest of health) that you:

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P2 respiratory protection.
- Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.
- Keep the muffle under water whilst devesting or
- Let the investment soak in water before devesting
- Liquid is alkaline!
- Liquid spilled on the bare skin should be washed off immediately using clean water and detergent!
- In case of contact with the eyes, rinse out **immediately** with water and call the nearest Doctor or Medical Centre!
- NEVER invest without metal casting ring!
- Never mix Deguvest liquid with other liquids
- Never mix Deguvest powder with other powders
- Keep liquid protected from frost

Deguvest SR should not be allowed to come into contact with *gypsum* or *materials containing gypsum*, as this could lead to *sulphur imposed damage* of certain metals with a high melting point! Therefore, keep your mixing utensils for **Deguvest SR** separate from others!

2 Technical Information

Deguvest SR

Mixing ratio (powder : liquid)	100 g : 22 ml
Mixing time (under vacuum)	60 s
Effective handling time	4 min – 6 min
Total expansion (linear) = setting expansion + thermal expansion	1,2 % – 2,4 %
Pressure resistance (depending on liquid concentration)	4 MPa – 8 MPa

3 Investment Storage

Most of you are probably aware of the problems which can occur when investment is subjected to variations in environment temperature (e.g. summer and winter) and the effect this can have on your castings. That is why it is very important to find the right place to store your supply of investment.

Deguvest SR should be stored in a cool and dry place, between 20 °C – 23 °C (71,6 °F – 73,4 °F).

If you have a place in your laboratory where these conditions can be found throughout the year, then that's the place to store your Investment! Under these conditions, we guarantee a shelf life of 24 months for both liquid and powder.

If it is not possible for you to store your investment at this temperature, then you should at least make sure that the work bench where investing is done fits these requirements (use a temperature control cupboard or similar).

Avoid any nasty surprises by observing the following:

- Avoid using any investment that has just been brought out of an especially warm or cold place. (e.g. storage rooms, silos, cupboards etc.) Fill up stocks at your workbench before your supply runs out. Give the investment time to adjust to the environment!
- Avoid using investment that has just been delivered to the laboratory. Again, give the investment time to acquire the right temperature.
- Sporadic temperature checks at the workbench help maintain constant temperatures.
- If space allows, stock up on investment (especially LIQUID) before the cold winter months arrive. This way, you forego the risk of exposing the liquid to frost during delivery.

4 Investment Instructions

Now you are ready and eager to get started and your first copings invested? Well, as a precaution, we would like to take you through a *few steps* which we know every Technician is aware of but, nevertheless, we feel worth mentioning.

4.1 Preparation

Our motto here is *"Shining utensils = shining results!"*

*Tip: Use a new flask right from the start. **Deguvest SR** helps keep mixing flask, spatula and environment exceptionally clean. You will notice also, **Deguvest SR** is not as hard on your hands as normal investments!*

- Take a look at your mixing flask (and the mixing arm itself)! Are they both clean? No residue or dried investment anywhere? Extremely old flasks should be replaced, as they can soak up between 2 and 3 ml of liquid when they are dry!

- Use a clean, dry spatula.
- If you use a wetting agent (e.g. DeguDent Waxit®), then ensure the surface of your copings are free from any residue before investing.
- Casting rings should be clean and lightly coated with vaseline.
- Check the position of your copings within the casting ring (height, distance from centre, sprues etc.).
- Use a slightly damp ring liner or alternatively fleece liner for speed heating. Fit the edge of the liner so that it is flush with the top rim of the casting ring. Secure with wax if necessary.

4.2 A few words about measuring out liquid

It is of utmost importance that **Deguvest SR** only be mixed using **Deguvest Liquid!** The total expansion of the investment can be varied by the adding distilled water. In this way, the specific shrinkage of the metal used (which occurs during cooling) can be compensated. The higher the liquid concentration, the more total expansion of the investment.

Table 2: Mixing ratios

Ring size	1 x ring	3 x ring	6 x ring	9 x ring
Powder	60 g	150 g	300 g	450 g
Liquid	13,2 ml	33 ml	66 ml	99 ml

We have prepared a chart to help you in measuring out the liquid (table 3). We have included examples of metals and the recommended dosage for each type in percent (%). The standard procedure at this point would be to work out the dosage for each individual investment in millilitres, measuring out your liquid in a small measuring cylinder. This method, however, does have its dangers:

- *It is not easy to measure out small amounts exactly.*
- *When converting from percent to millilitres it can occur that you reach figures with a decimal point (e.g. 28,05 ml). Such a number is impossible to measure out in a cylinder.*
- *Dried out liquid invariably leaves a crust inside the measuring cylinder. If cylinders are not regularly replaced, they can influence measurements.*

We therefore recommend a different approach to measuring out your liquid. Use a large calibrated 1 litre glass cylinder to measure out the percentage you require for your daily investments. For example, if you use a 75% concentration for your gold framework, then you can measure out 750 ml liquid to 25 ml distilled water and then fill this mixture into a plastic bottle. All you need to do then is measure out the required amount in ml for your investment. Using the large cylinder you can measure out more exact and you also have the same concentration over a long period (the life of your 1 litre bottle)

4.3 Mixing powder and liquid

Mixing Ratio (powder : liquid) = 100 g : 22 ml

- Once measured out, pour the liquid into the clean mixing flask.
Important! Older flasks with a rough inner lining can soak up between 2-3ml of liquid when they are dry! Before using such a flask, you should always rinse out with water and dry off lightly with a towel.
- Weigh out the required powder or alternatively use the prepacked sachets and add to the liquid. (Be careful not to spill any powder as this would obviously change the mixing ratio and thus your results will not be as expected!)
- Mix the two components manually for **30 seconds**, making sure that no dry powder remains on the flask wall.
Close the mixing flask. If your mixer allows (except using **DeguDent Multivac® compact**, see below) then the flask should be placed under vacuum for a further **30 seconds**. Always make sure that you have the **highest level of vacuum** possible. This determines the quality of your results!
Mix for **60 seconds** under vacuum in the mixer (Multivac compact).

Be aware of your rpm!

The revolutions per minute of your unit's mixing arm has a direct effect on the setting time / setting expansion of the investment. The technical information given in these instructions are all based upon tests using a constant 360 rpm. Should it be that your mixing unit varies from this, then it may cause slight variations in casting results.

If working with the **Multivac compact**, please use program no.1 (Deguvest CF). The premixing of powder and liquid is not necessary. Alternatively, use the standard program and proceed as described above.

4.4 Investment

Using a **light** vibration, **Deguvest SR** flows easily into the casting ring and copings. Avoid trapping air around difficult areas on the coping. After the wax coping is covered, turn off the vibrator and continue filling the casting ring up to the rim.

Due to the smooth consistency and a handling time of 4 to 6 minutes you will have no problem investing even the most complicated structures or copings.

If you use a **Multivac compact** unit, then you should call up program no. 1 (Deguvest CF) or the standard setup and follow the instructions for this mixing unit.

5 Setting times

If you wish to use the speed or quick cast method, then the casting ring must be placed in the furnace at casting temperature **exactly 15 minutes after investing!** We recommend setting a stop watch or timer to 15 minutes and pressing start the moment you start mixing. When the 15 minutes have elapsed, you can place the ring straight in the furnace. If you wish to heat your casting ring in the conventional way, then we recommend you use the following setting times before placing in the cold furnace:

Size 1 ring = 20 min

Size 3 ring = 30 min

Size 6 ring = 40 min

Size 9 ring = 50 min

6 Preheating and casting

6.1 Conventional Heating

After the investment has been allowed to set, place the casting ring in the COLD furnace. (Cold means below 100°C/212°F). The heating rate should be 7°C/min (44,6°F/min) and the firing program should include the following steps (depending on the size of the ring):

- 45 to 60 minutes at 270°C/518°F (Cristobalite transformation)
- 30 to 60 minutes at 560°C/1040°F (Quartz transformation)
- holding time at final temperature:

Size 1 ring = 20 min

Size 3 ring = 30 min

Size 6 ring = 40 min

Size 9 ring = 50 min

6.2 Quick Cast Method

For **speed or quick cast** method, please observe the following:

- Always use a metal casting ring!
- Do not exceed a size 6 casting ring!
- Always use ceramic liner (Deguvest fleece)!

We recommend the muffle should not be placed flat on the furnace floor, rather it should be tilted slightly allowing the wax to flow freely from the casting ring. A good method of achieving the required tilt is to place the casting ring on a piece of honeycombed firing carrier as used in ceramics!

CAUTION: DURING THE FIRST FEW MINUTES OF BURNOUT IN THE FURNACE TOXIC GASES AND FLAMES ARE PRODUCED. AVOID OPENING THE FURNACE DOOR UNTIL THE HOLDING TIME HAS ELAPSED!

Holding times at firing temperature must always be followed!

Size 1 ring = 20 min

Size 3 ring = 30 min

Size 6 ring = 40 min

Please take into account that certain parameters outside of those mentioned here, which include material (wax, composites, etc.), geometrical forms and size of copings, casting ring size and also type of alloy being used, do have an effect on the results achieved.

We have arranged all this information for you in a chart which you will find as enclosure with these instructions (page 24).

7 Casting

Please observe the casting instructions given for the alloys which you use, as well as the instructions for the use of casting units.

8 Devesting

Caution! Deguvest SR loses its dust reducing attributes after it has been heated in the furnace!

Allow the casting ring to cool at room temperature until hand warm. We recommend soaking the ring in water approximately 15 minutes before devesting. Not only is the casted framework easier to free from investment, it also reduces dust development whilst working. Use devesting pliers to free the framework from most of the investment. Following this, the framework can be sand blasted with aluminium oxide (250µm) or alternatively cleaned by placing the object for 15 minutes in pickling agent, then for a further 15 minutes in an ultrasonic cleaner. Any remaining investment can then be removed using a toothbrush and water. (The surface of the object profits from this method).

9 Recommendations for Non-precious alloys

Deguvest SR investment is stable up to 1000 °C and compatible with all widespread non-precious alloys (except for titanium) – and is also suitable for speed casting. Please observe the alloy manufacturers' recommended preheating temperatures.

We advise using a liquid concentration of 80–90% for single copings and 85–95% for bridge restorations. Because of the sheer number of dental alloys available and their varying characteristics, we advise making reference castings before attempting complex restorations.

We wish you much success with **Deguvest SR** investment. If you have any questions as regards handling or technical details, please do not hesitate to call your local DeguDent representative, or visit our website at www.degudent.com.

Your **DeguDent** Team

Investing with Deguvest SR

Quick reference guide

- 1 Mixing ratio: 100 g : 22 ml
- 2 Prepare utensils, cleaning (mixing flask, mixing arm, spatula).
- 3 Clean and grease sprue former, casting ring. Prepare ring liner or fleece.
- 4 Attach sprues as normal.
- 5 If using wetting agent, blow dry lightly. No residue!
- 6 Measure out liquid required (table 1, chapter 4.2) Temperature?
- 8 Weigh out the powder. Temperature?
- 9 Fill the flask with both components. If Quick Casting, start your timer now!
- 10 Pre-mix per hand for 30 seconds.
- 11 Close flask. 30 seconds vacuum only!
- 12 Start mixing unit. (60 seconds).
- 13 Investment begins.
- 15 If Quick Casting : remember to place the ring in the furnace straight after the 15 minutes have elapsed!

Tip: For easier reference, just tear this page out of the booklet and pin up at your place of work.

Table 3: Mixing ratios of various liquid concentrations

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Tabelle 4: Typical concentrations of Deguvest Liquid for investing objects in Deguvest SR

		High gold content alloy	Reduced gold content	Pd-based alloy	Notes
Occlusal Inlay	wax	55 % - 65 %			
MOD Inlay	wax	60 % - 70 %			
Anterior crown	wax Adapta	60 % - 70 %	65 % - 75 %	70 % - 80 %	
Posterior crown	wax Adapta	60 % - 70 %	65 % - 75 %	70 % - 80 %	
3-4 unit bridge	wax Adapta	60 % - 70 %	70 % - 80 %	70 % - 80 %	
bridge above 6 units	wax Adapta	45 % or less	45 % or less	45 % or less	
0° Telescope second.*	Adapta Pattern	60 % - 70 % 45 % - 55 %			
6° secondary*	Adapta Pattern	50 % - 60 % 40 % - 50 %			
half circular milled attachment with shoulder/pin	Pattern	45 % - 55 %			

* We do not recommend the rapid heating method for preheating complex restorations.

Please take into account, that a large number of parameters (e.g. temperature, modelling materials, mixing unit etc.) have a direct effect on the results and that the above recommendations may therefore vary individually.

Tabelle 5: Heat treatment: Conventional Heating

Size of the casting ring	1 x	3 x	6 x	9 x
Setting time	20 min	30 min	40 min	50 min
Furnace start temperature	<100° C	<100° C	<100° C	<100° C
Heating rate	7°C / min	7°C / min	7°C / min	7°C / min
1. Holding stage	270 °C	270 °C	270 °C	270 °C
Holding time	20 min	30 min	40 min	50 min
Heating rate	7°C / min	7°C / min	7°C / min	7°C / min
2. Holding stage	580 °C	580 °C	580 °C	580 °C
Holding time	20 min	30 min	40 min	50 min
Heating rate	7°C / min	7°C / min	7°C / min	7°C / min
3. Holding stage	final temperature	final temperature	final temperature	final temperature
Holding time.	20 min	30 min	50 min	60 min
Casting				

Tabelle 6: Heat treatment: Rapid Heating

Size of the casting ring	1 x	3 x	6 x	9 x
Setting time	15 min	15 min	15 min	not recommended
Preheated furnace temperature	final temperature	final temperature	final temperature	not recommended
Holding time at final temperature	30 min	40 min	50 min	not recommended
Casting				not recommended

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DeguDent GmbH