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# Optivest<sup>®</sup> avanti

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso  
Инструкция по  
применению

**DeguDent**

*A Dentsply International Company*

# Optivest® avanti

<b>D</b>	Gebrauchsanweisung . . . . .	3
<b>GB</b>	Instructions for use . . . . .	7
<b>F</b>	Mode d'emploi . . . . .	11
<b>I</b>	Istruzioni per l'uso . . . . .	15
<b>E</b>	Instrucciones para el uso . . . . .	19
<b>RUS</b>	Инструкция по применению . . . . .	23

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Optivest avanti is a phosphate bonded special investment compound that can be heated quickly and is suitable for the manufacturing of duplicate models and casting moulds for the CoCr model casting technique. Optivest avanti does not contain any plaster or graphite. After drying the casting mould can be placed directly into the preheating furnace that has been heated up to 900 °C. Optivest avanti is mixed with Special liquid (REF 2540 0114).

## 1 Technical data

### 1.1 Characteristic material data according to DIN 13 919

(with Special Liquid in a concentration of 50 %)

Setting time according to Vicat:	6 – 8 min
Setting expansion:	< 1 %
Thermal expansion:	ca. 0.8 %
Compressive strength:	> 15 MPa

### 1.2 Processing recommendations

Material temperature:	21 °C
Mixing ratio:	100 g : 20 ml
Mixing time under vacuum:	60 s
Working time range:	5 – 7 min

## 2 Important information

Prior to use observe the safety instructions provided under 6 and 6.1. Optivest must be stored at a dry and cool place. If stored properly, the storage time of the investment material amounts to 18 months. The Special Liquid must be stored at a room temperature between 5 °C and max. 30 °C.

### **3 Expansion control**

The total expansion, i.e. the sum of the setting expansion and the thermal expansion can be changed by diluting the Special Liquid with distilled water. Depending on the alloy that is used, the concentration of the Special Liquid is between 45 % and 55 %.

## **4 Manufacturing of duplicate models**

### **4.1 Mixing of the investment material**

Premix the powder and the liquid with a spatula in a clean, freshly rinsed bowl of the Multivac® evacuating unit. Then let the investment material under vacuum for a short time and mix under vacuum for 60 sec. Please use program number 3 for mixing with the Multivac compact. For the manufacturing of investment models we recommend a concentration of the mixing liquid of 45 – 55 %.

### **4.2 Pouring up of the silicone duplicating mould**

Silicone duplicating moulds can be poured up more easily if they have previously been wetted with Waxit®. The moulds must be blown dry before pouring up with investment material (do not use compressed air!). Let the investment material flow into the duplicating mould while vibrating slightly and avoid the formation of bubbles during this process. After the duplicating mould has been filled, switch off the vibrator immediately.

### **4.3 Pouring up of agar-agar duplicating moulds**

If duplicating gel (Geloform) is used in combination with Optivest avanti, it is recommended to soak the master model in luke-warm water for approx. 10 min prior to duplicating.

### **4.4 Setting**

The investment model should set at least 30 min; then it is carefully removed from the duplicating mould. In order to prepare undercuts which result in an improved connection of model /coat, the investment model should be provided with oblique edges using the trimmer (dry!).

### **4.5 Drying**

Dry the model for approx. 30 min in a drying cabinet which has been preheated to 150 °C. Then apply adhesion agent to improve adhesion of the wax pattern.

## 4.6 Hardening

Optivest avanti investment models which have been made with silicone duplicating moulds do not require hardening. If the investment model has been produced with duplicating gel, hardening with cold hardener is recommended. Please observe the respective information for use of the corresponding product.

# 5 Manufacturing of the casting mould (coating)

## 5.1 Reduction of tension

Apply or spray a thin layer of Waxit wetting agent to the wax object and then blow it dry carefully (do not use compressed air!).

## 5.2 Fine investing

Due to the fine-grained composition of Optivest avanti fine investing is not required.

## 5.3 Mixing the investment compound

Dose Optivest avanti powder and Special liquid and mix as described under 1.1. We recommend to select a concentration of the Special liquid of 45 – 55 %. The mixing ratio of water : liquid should be 100 g : 20 ml.

## 5.4 Investing

Let the mixed investment material flow into the casting mould while vibrating slightly; the formation of bubbles must be avoided. After the casting mould has been filled, switch off the vibrator immediately.

## 5.5 Setting

Let the casting mould set for approx. 15 min. After the investment compound has set, the upper surface at the casting ring should be trimmed slightly to remove the glaze layer so that the escape of gas is not inhibited.

## 5.6 Preheating

After setting, the casting mould can be immediately placed into the preheating furnace that has been heated up to 900 °C. Then heat up resp. cool down to the casting temperature recommended by the manufacturer of the alloy. At the preheating temperature of 900 °C the investment material performs an additional crystallization process which ensures extremely high temperature resistance.

## 5.7 Casting

Cast according to the instructions provided by the manufacturer of the alloy.

## 6 Safety information

Mixing with other liquids or materials during this process health hazards as well as impairments to the cast quality cannot be excluded.

All dental investment material contain parts of crystalline quartz and /or cristobalite.

These substances can adversely affect health if they are breathed in. Can cause silicosis. Always use a sufficiently powerful extractor or wear a mask.

Carry out the “coarse” removal of the casting under water or soak the investment material mold.

### 6.1 Information for the processing of mixing liquids for dental investment materials

Following contact with skin:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

The Safety Data Sheet includes further details.

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