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Silber Star

Instructions for use

SilberStar

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SilberStar is a plaster-free, phosphate bound special investment compound for the production of duplicate models and casting moulds for the CoCr model casting technique. SilberStar is particularly suitable when used in combination with silicone duplicating materials. The fine-grained SilberStar powder produces extremely smooth surfaces and thus facilitates finishing considerably. The matched relationship between expansion and strength results in restorations with high precision of fit. SilberStar is mixed with a special mixing liquid with properties that match the properties of SilberStar powder perfectly.

1 Technical data

Mixing ratio powder : liquid	100 g : 15 ml
Total expansion	0.7 – 1.8 %
Mixing time under vacuum	60 sec
Processing time	3 – 4 min
Compressive strength depending on the concentration of the mixing liquid	15 – 20 MPa
Recommended temperature for processing and storage	16 – 20 °C

2 Important information

Prior to use observe the safety instructions provided under 6 and 6.1 SilberStar must be stored at a dry and cool place. If stored properly, the storage capacity of the investment compound including the liquid amounts to 18 months. SilberStar mixing liquid must be stored at a room temperature between 5 °C and max. 30 °C.

3 Expansion control

The total expansion, i.e. the sum of the setting expansion and the thermal expansion can be changed by diluting the mixing liquid with distilled water.

Concentration of the mixing liquid (%)	Setting expansion (%)	Thermal expansion (%)	Total expansion (%)
0	0.1	0.6	0.7
50	0.2	0.9	1.1
75	0.4	1.1	1.5
100	0.6	1.2	1.8

4 Production of duplicate models

4.1 Mixing of the investment compound

Premix powder and liquid in a clean mixing bowl of the evacuating unit. Then place the investment under vacuum for a short time and mix under vacuum for 60 sec. For the production of investment material models we recommend a concentration of the mixing liquid of 75 %.

4.2 Pouring up of the silicone duplicating mould

Silicone duplicating moulds can be poured up more easily if they have previously been wetted with Waxit®. The moulds must be blown dry before pouring in the investment. Let the investment material flow into the duplicating mould while vibrating lightly; avoid the formation of bubbles during this process. After the duplicating mould has been filled, switch off the vibrator immediately.

4.3 Pouring up of agar-agar duplicating moulds

If duplicating gel is used in combination with SilberStar, it is recommended to place the master model in luke-warm water for approx. 10 min prior to duplicating. Recommended concentration for the investment material model: 50 %.

4.4 Setting

The investment material should set for at least 30 min; then it is carefully removed from the duplicating mould. In order to prepare undercuts which result in an improved connection of model/coat, the investment material model should be provided with oblique edges using the trimmer.

4.5 Drying

Dry the model in a drying cabinet which has been preheated to approximately 80 °C for approx. 40 min. Then apply adhesion agent to improve adhesion of the wax pattern.

4.6 Hardening

SilberStar models which have been produced in conjunction with silicone duplicating moulds do not have to be hardened. If the investment compound model has been produced with duplicating gel, hardening with cold hardener is recommended. Please observe the respective instructions for use of the corresponding product.

5 Preparation of the casting mould (lining)

5.1 Reduction of tension

Apply or spray a thin layer of wetting agent to the wax object and then blow it dry carefully (do not use compressed air!).

5.2 Fine investing

Due to the fine grain size of the SilberStar powder and the compatibility with all modelling materials, fine investing is not required.

5.3 Mixing the investment compound

Dose and mix SilberStar powder and mixing liquid as described in 4.1. We recommend to select a concentration of the mixing liquid of 60%. The mixing ratio of water : liquid should be 100 g : 15 ml.

5.4 Investing

Let the mixed investment compound flow into the casting mould under light vibration; the formation of bubbles must be avoided. After the casting mould has been filled, switch off the vibrator immediately.

5.5 Setting

Let the casting mould set for approx. 60 min. After the investment compound has set, the upper surface at the casting ring should be trimmed slightly in order to remove the glaze layer to ease the escape of gases.

5.6 Wax elimination

Place the casting mould into a furnace which has been preheated to 300 °C and eliminate the wax for approx. 45–60 min. Place the casting ring into the furnace with the casting funnel facing downward.

5.7 Preheating

After the wax has been eliminated, bring the casting mould to the final temperature of 1000–1050 °C recommended by the manufacturer. The heat rate should not exceed 8 °C/min. The final temperature should be kept for 60 min; if a larger number of casting rings is used, the time must be extended correspondingly.

5.8 Casting

Cast according to the processing instructions provided by the manufacturer of the alloy.

6 Safety information for the processing of dental investment compounds

Mixing with other liquids or compounds must not be performed since health hazards as well as negative influences on the quality of the cast objects cannot be excluded.

All dental investment compounds contain proportions of crystalline quartz and/or cristobalite. These substances can adversely affect health if they are inhaled. Always use a sufficiently powerful extractor or wear a mask.

“Coarse devesting” after casting is performed under water.

If heated over 200 °C: formation of ammonia. Ammonia irritates the respiratory organs.

6.1 Information for the processing of mixing liquids for dental investment compounds

Following contact with skin:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

Further information is contained in the safety data sheet.

U3 (S.11)

Vakat !

www.degudent.com

DeguDent GmbH
P.O. Box 13 64 - 63403 Hanau
GERMANY
www.degudent.com